

Pipettes straight



Pipettes, straight, serological, glass, graduated, non-sterile



02

Non-sterile packed in plastic sleeves of 10 within the standard pack. Accuracy is within $\pm 2\%$ of volume, which is adjusted for blowout. 10mL size has polypropylene top stem of reduced diameter, this enables finer control to be achieved when hand pipetting and allows superior plugging. Black graduations. PMK-100-020A/030U supplied with PVT test and adapter.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Length, mm	Jet	Pack qty
PMK-100-020A	D804	1	0.01	230	Short	500
PMK-100-030U	D805	1	0.10	230	Short	500
PMK-100-040R	D802	1	0.01	270	Drawn	500
PMK-100-050X	D806	2	0.02	260	Short	500
PMK-100-060L	D803	2	0.02	270	Drawn	500
PMK-100-070Y	D807	5	0.05	310	Drawn	500
PMK-100-080F	D808	10	0.10	310	Drawn	250



Pipettes, straight, serological, glass, graduated

02



Sterile in inner sleeves of ten; 5mL and 10mL in inner sleeves of five.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Length, mm	Jet	Pack qty
PMK-200-040J	D804/PS	1	0.01	230	Short	500
PMK-200-041H	D805/PS	1	0.10	230	Short	500
PMK-200-042F	D802/PS	1	0.01	270	Drawn	500
PMK-200-050G	D806/PS	2	0.02	260	Short	500
PMK-200-052C	D803/PS	2	0.02	270	Drawn	500
PMK-200-070A	D807/PS	5	0.05	310	Drawn	500
PMK-200-080U	D808/PS	10	0.10	310	Drawn	250



Pipettes, straight, serological, with magnifier stripe



Suitable for tissue culture work.

A black magnifier stripe runs lengthwise along the pipette making the meniscus obvious and reducing reading errors. The pipettes have excellent optical clarity with high contrast black markings making it easier to read volumes. The 5mL, 10mL, 25mL and 50mL sizes have both ascending and descending graduations. Tip design on the 25mL and 50mL sizes increases the effect of surface tension to minimise dripping. The individually wrapped pipettes come in a clear, all plastic wrap and the multipacks come in a recyclable LDPE bag.

Individually wrapped

Catalogue No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
FB55481	1	0.01	Yellow	100	1,000
FB55482	2	0.01	Green	100	500
FB55483	5	0.1	Blue	50	200
FB55484	10	0.1	Orange	50	200
FB55485	25	0.2	Red	50	200
FB55486	50	0.5	Purple	25	100

Multipack

Catalogue No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
FB55488	1	0.01	Yellow	25	1,000
FB55489	2	0.01	Green	25	500
FB55490	5	0.1	Blue	25	500
FB55491	10	0.1	Orange	25	500
FB55492	25	0.2	Red	20	200
FB55493	50	0.5	Purple	20	100