

Pipettes, straight, serological, polystyrene, paper/plastic wrap, Stripette®

CORNING

66



Ideal for sterile tissue culture applications, especially while wearing latex gloves as the wrap reduces static cling.

A strip creates a thermometer-like effect at the meniscus, making the volume easy to read. Colour coding on wrapping allows easy size identification. Certified non-pyrogenic. Individually packed.

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-670-011G	4485	1	Yellow	50	1,000
TKV-670-031A	4486	2	Green	50	1,000
TKV-670-051R	4487	5	Blue	50	200
TKV-670-071L	4488	10	Orange	50	200
TKV-670-091F	4489	25	Red	25	200
TKV-670-111C	4490	50	Purple	25	100
TKV-670-121W	4491	100	Aqua	10	100



Pipettes, straight, serological, polystyrene, plastic wrap, bulk pack, Stripette®

CORNING


66



Bulk pack pipettes are ideal for large scale sterile and non-sterile liquid handling

A strip which creates a thermometer-like effect at the meniscus, making the volume easy to read. Certified non-pyrogenic.

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-673-011L	4010	1	Yellow	50	1,000
TKV-673-031F	4020	2	Green	50	1,000
TKV-673-051W	4050	5	Blue	50	500
TKV-673-071Q	4100	10	Orange	50	500
TKV-673-081N	4250	25	Red	50	200
TKV-673-091K	4500	50	Purple	50	100



Pipettes, straight, serological, polystyrene, plastic wrap, Stripette®

CORNING

66



Ideal for sterile applications such as tissue culture.

Feature a unique vertical strip which creates a thermometer-like effect at the meniscus, making volume easy to read. 'Banana-peel' plastic wrap provides easy, aseptic opening and pipette removal. Individual, clear, plastic wrap allows full view of pipette for easy size identification and reduces storage space requirements. Certified non-pyrogenic. TKV-671-090B has a wide tip for handling viscous fluids.

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-671-010C	4011	1	Yellow	100	1,000
TKV-671-020W	4012	1	Yellow	100	200
TKV-671-030T	4021	2	Green	100	1,000
TKV-671-040Q	4051	5	Blue	50	200
TKV-671-050N	4101	10	Orange	50	200
TKV-671-060K	4492	10	Orange	50	200
TKV-671-070H	4251	25	Red	50	200
TKV-671-080E	4501	50	Purple	25	100
TKV-671-090B	4484	100	Aqua	10	100

