

Pipettes, LDPE, Pasteur

41

Fisherbrand®
S

Single piece polyethylene moulded pipette.
Maximum bulb capacity, 3.5mL. Approx. length 153mm.

Catalogue No	Description	Graduations, mL	Inner pack qty	Pack qty
FB55343	1mL non-sterile	0.25	-	500
FB55345	3mL non-sterile	0.5	-	500
FB55340	1mL sterile	0.25	1	500
FB55348	3mL sterile	0.5	1	500
FB55341	1mL sterile	0.25	10	500
FB55349	3mL sterile	0.5	10	500
FB55358	Fine point, non-sterile	-	-	500
FB55359	Thin stem, non-sterile	-	-	500



Pipettes, polyethylene, transfer, Pastette®

02

S

Precise, accurate and unbreakable. Does not require a separate bulb. Non-toxic to cells and microorganisms, with a low affinity for proteins.



Graduated Pastette®

Catalogue No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
PMK-500-010R	LW4212	Small bulb graduated Pastette®	138	1.5	25	40	0.25	500
PMK-500-020X	LW4040	1mL	150	3.5	25	40	0.25	500
PMK-500-030L	LW4005	1mL sterile 1's	150	3.5	25	40	0.25	500
PMK-500-040Y	LW4006	1mL sterile 10's	150	3.5	25	40	0.25	500
PMK-500-050F	LW4010	1mL sterile 20's	150	3.5	25	40	0.25	500
PMK-500-060C	LW4274	Medi-Pastette®	156	2.0	25	40	0.50	500
PMK-500-070W	LW4111	3mL	156	3.5	22	45	0.50	500
PMK-500-080T	LW4112	3mL sterile 1's	156	3.5	22	45	0.50	500
PMK-500-090Q	LW4113	3mL sterile 10's	156	3.5	22	45	0.50	500
PMK-500-100Q	LW4114	3mL sterile 20's	156	3.5	22	45	0.50	500
PMK-500-110N	LW4282	Mini-Pastette®	114	1.0	23	43	0.1	500

Fine tip Pastette®

Catalogue No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
PMK-510-010K	LW4060	Fine tip	144	3.5	50	20	-	500
PMK-510-020H	LW4061	Fine tip sterile 1's	144	3.5	50	20	-	400
PMK-510-030E	LW4070	Fine tip sterile 20's	144	3.5	50	20	-	500
PMK-510-040B	LW4635	Extended	153	3.5	50	20	-	500
PMK-510-050V	LW4234	Maxi-Pastette®	150	6.0	45	22	-	400
PMK-510-060S	LW4235	Maxi-Pastette® sterile 1's	150	6.0	45	22	-	400
PMK-510-070P	LW4236	Maxi-Pastette® sterile 20's	150	6.0	45	22	-	400
PMK-510-080M	LW4231	Extended Mini-Pastette®	104	1.0	50	20	-	400
PMK-510-090J	LW4233	Extended	153	3.0	50	20	-	250
PMK-510-100J	LW4232	Extended sterile 20's	153	3.0	50	20	-	500

Micro Pastette®

Catalogue No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
PMK-510-110G	LW4237	Extended fine tip	153	3.0	50	20	-	500
PMK-520-020A	LW4242	Short stem	83	3.5	25	40	-	750
PMK-520-025N	LW4858	Short stem non sterile	60	-	18	56	-	500
PMK-520-030U	LW4020	Standard	150	3.5	25	40	-	500
PMK-520-040R	LW4252	Large bulb	150	7.7	22	45	-	250
PMK-520-050X	LW8241	Short-stem ultra	58	1.0	25	40	-	500
PMK-520-060L	LW4243	Fine tip ultra	50	1.0	45	22	-	500

Small Pastette®

Catalogue No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
PMK-530-010T	LW4273	3 inch	87	1.0	25	40	-	500
PMK-530-020Q	LW4295	Paddle	124	0.2	20	50	-	500
PMK-530-030N	LW4293	Fine tip paddle	130	0.2	33	30	-	500

entry continued