

Pipette fillers, Pi-pump

03

An almost indestructible acid resistant pipette controller, easily dismantled for cleaning. Suitable for filling and dispensing from all standard pipetts from micro to 25mL.

Catalogue No	Colour	For pipette size, mL
PMR-336-005X	Yellow	0.2
PMR-336-010V	Blue	2
PMR-336-030P	Green	10
PMR-336-050J	Red	25



Pipette fillers

41

Fisherbrand®

Moulded rubber bulb, approx. 60mL capacity, with three glass ball valves which can be controlled by finger pressure. No metal parts to corrode and can be used for all liquids except those which attack rubber. Can be used on bulb or cylindrical pipettes from 2mL to 50mL.
Note: Do not allow liquids to come in contact with the valves or inside of the bulbs - this may impare performance.

Catalogue No	Colour
FB51420	Black
FB51421	Blue
FB51422	Red
FB51423	Yellow



Pipette fillers, PipetBoy Acu

04

IBS INTEGRA BIOSCIENCES
Designed for the Scientist

S

Cable free, pipette controller.
Extremely lightweight, compact and easy to use. Works on a rechargeable 9V NiMH battery which is easily inserted at the top of the handle. One charge is sufficient to fill and purge at least 500 x 50mL pipettes. Accepts standard pipettes ranging from 0.1mL to 100mL. Comes in a range of five transparent colours and a traditional light grey to aid identification. With battery, transformer for 230V/50Hz, sterile filter, wall mount and a Swiss army knife.

Catalogue No	Alt. No	Colour
PMR-100-010C	155000	Light grey
PMR-100-030T	155015	Green
PMR-100-050N	155016	Red
PMR-100-070H	155017	Blue
PMR-100-090B	155018	Clear
PMR-100-095X	155019	Purple

Accessories

Catalogue No	Alt. No	Description
PMR-100-510F	153005	Rechargeable battery 9V NiMH
PMR-100-520C	151005	Handle mount
PMR-100-530W	151020	Pipettor, rubber, 0.1mL to 50mL
PMR-100-115L	153015	Filter, sterile, for PipetBoy

