




Bulk pack



Racked open



Racked sterile



Reload

## Pipettor tips, universal fit

### SureOne Rack

- Flush Fit Inserts sit flush with rack, allowing easy multichannel pipetting
- Notched clasp designed to prop lid up for autoclave venting
- Sure Stack Lid rim prevents racks from sliding around
- Safe Lock Latch allows for one-handed opening and closing ability and won't pop open when dropped
- Sure Snap design securely holds insert in place, preventing it from being lifted out while pipetting
- Flip-Lid can be either a hinged or lift off lid; hinged format does not prevent use of multichannel pipettor

### SureOne Insert

- Contact Clip prevents unnecessary contact with the tips while clipping in inserts
- Lot Numbered Inserts for traceability include tip volume and Fisher Catalogue No
- Indented Insert Holes inhibit tip movement, stabilising them for multichannel work
- Rigid Inserts are all 96 well multichannel format offering a robust non-pliable design providing secure tip seal

### SureOne Reload Packaging

- Dust Flap prevents debris from contaminating the product
- Canyon Cut Package allows for easy tip insert removal

A comprehensive range of non-filter and aerosol barrier pipettor tips with the environmentally friendly reload system equalling less waste.

Pipettor tips have a universal fit and are compatible with most popular brands of pipettor.

### Standard Sterility Claims

Non sterile products - certified to be free of RNase/DNase and DNA

Sterile products - e-beam sterile products are certified to be free of RNase/DNase, Pyrogen, Bioburden, PCR\* inhibitors, Endotoxin free.

### Filter tips Racked sterile

Catalogue No	Description	Volume, $\mu\text{L}$	Colour	Inner pack qty	Pack qty
1 FB78098	Micropoint tip, graduated at 2 $\mu\text{L}$	0.1 to 10	Clear	96	960
FB78100	Micropoint tip, graduated at 2.5 $\mu\text{L}$ , extended length	0.1 to 10	Clear	96	960
FB78102	Micropoint tip Eppendorf style	0.1 to 10	Clear	96	960
FB78104	Micropoint tip Eppendorf style	2 to 20	Clear	96	960
2 FB78106	Bevelled tip	2 to 20	Clear	96	960
FB78108	Bevelled tip	10 to 100	Clear	96	960
FB78110	Bevelled tip, graduated at 10, 50, 100 $\mu\text{L}$	20 to 200	Clear	96	960
FB78112	Micropoint tip, graduated at 100, 200, 500, 1,000 $\mu\text{L}$	100 to 1,000	Clear	96	960

### Standard tips

Micropoint tip, graduated at 2 $\mu\text{L}$ . 0.1 $\mu\text{L}$  to 10 $\mu\text{L}$

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
FB78000	Bulk bag	Clear	1	1,000
FB78002	Racked	Clear	96	960
FB78004	Racked sterile	Clear	96	960
FB78006	Reload system pack	Clear	96	960

Extended length Micropoint tip, graduated at 2.5 $\mu\text{L}$ . 0.1 $\mu\text{L}$  to 10 $\mu\text{L}$

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
FB78008	Bulk bag	Clear	1	1,000
FB78010	Racked	Clear	96	960
FB78012	Racked sterile	Clear	96	960

Micropoint tip Eppendorf style. 0.1 $\mu\text{L}$  to 20 $\mu\text{L}$

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
FB78014	Bulk bag	Clear	1	1,000
FB78016	Racked	Clear	96	960
FB78018	Racked sterile	Clear	96	960



\*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche