

Pipettors, fixed volume, Thermo Scientific Finnpiquette® F2

04



- Advanced volume gearing mechanism to improve accuracy and precision
- Fully autoclavable
- Double action pipetting button with rotating upper part independent of volume setting

- Large ergovision display
- Light pipetting action
- 'Super Blow Out' mechanism and soft touch tip ejection



Catalogue No	Alt. No	Volume, µL	Finntip	Inaccuracy, %	Imprecision, cv%
PMP-025-005C	4652000	1	Flex 10, 10, 50	±4.0	4
PMP-025-010J	4652010	5	Flex 10, 10, 50	±1.4	1.4
PMP-025-015W	4652020	10	Flex 200, 250 Univ, 200 Ext	±0.9	0.8
PMP-025-020G	4652030	25	Flex 200, 250 Univ, 200 Ext	±0.6	0.5
PMP-025-025T	4652040	50	Flex 200, 250 Univ, 200 Ext	±0.6	0.4
PMP-025-030D	4652050	100	Flex 200, 250 Univ, 200 Ext	±0.4	0.3
PMP-025-035Q	4652060	250	Flex 1000, 1000, 1000 Ext	±0.4	0.3
PMP-025-040A	4652070	500	Flex 1000, 1000, 1000 Ext	±0.3	0.3
PMP-025-045N	4652080	1,000	Flex 1000, 1000, 1000 Ext	±0.3	0.3
PMP-025-050U	4652090	2,000	5mL	±0.3	0.2
PMP-025-055K	4652100	3,000	5mL	±0.3	0.2
PMP-025-060R	4652110	5,000	5mL	±0.3	0.2
PMP-025-065H	4652120	10,000	10mL	±0.3	0.2

Accessories

Catalogue No	Alt. No	Description
PMP-010-500G	9420400	Pipettor F-stand, 6 position

Pipettors, variable volume, Thermo Scientific Finnpiquette® F2

04



- Advanced volume gearing mechanism to improve accuracy and precision
- Fully autoclavable
- Double action pipetting button with rotating upper part independent of volume setting
- Large ergovision display
- Light pipetting action
- 'Super Blow Out' mechanism and soft touch tip ejection



Technical Specification - Specific

	PMP-020-005J	PMP-020-010Q	PMP-020-015G	PMP-020-020N	PMP-020-025D	PMP-020-030K
Volume, µL	2; 1; 0.2	5; 2.5; 0.5	10; 5; 1	10; 5; 1	20; 10; 2	20; 10; 2
Inaccuracy, %	± 2.5; ± 4.0; ± 12.0	± 1.5; ± 2.5; ± 6.0	± 1.0; ± 1.5; ± 2.5	± 1.0; ± 1.5; ± 3.5	± 1.0; ± 1.5; ± 3.0	± 1.0; ± 1.5; ± 3.0
Imprecision, cv%	2.00; 3.50; 10.00	1.00; 1.50; 5.00	0.50; 0.80; 2.00	0.80; 0.80; 3.00	0.40; 0.60; 2.50	0.40; 0.60; 2.50

	PMP-020-035A	PMP-020-040H	PMP-020-045U	PMP-020-050E	PMP-020-055R
Volume, µL	100; 50; 10	200; 100; 20	1,000; 500; 100	5,000; 2,500; 500	10,000; 5,000; 1,000
Inaccuracy, %	± 0.8; ± 1.2; ± 3.0	± 0.6; ± 1.0; ± 1.8	± 0.6; ± 0.8; ± 1.0	± 0.5; ± 0.7; ± 2.0	± 0.5; ± 0.8; ± 2.0
Imprecision, cv%	0.20; 0.40; 1.00	0.20; 0.40; 0.70	0.20; 0.30; 0.60	0.20; 0.30; 0.80	0.20; 0.30; 0.80

Catalogue No	Alt. No	Volume range, µL	Increments, µL	Colour code	Finntip
PMP-020-005J	4642010	0.2 to 2	0.002	Pink	Flex 10, 10, 50
PMP-020-010Q	4642020	0.5 to 5	0.01	Pink	Flex 10, 10, 50
PMP-020-015G	4642030	1 to 10	0.02	Pink	Flex 10, 10, 50
PMP-020-020N	4642040	1 to 10	0.02	Yellow	Flex 200, 250 Univ
PMP-020-025D	4642050	2 to 20	0.02	Turquoise	50
PMP-020-030K	4642060	2 to 20	0.02	Yellow	Flex 200, 250 Univ
PMP-020-035A	4642070	10 to 100	0.2	Yellow	Flex 200, 250 Univ, 300, 200 Ext
PMP-020-040H	4642080	20 to 200	0.2	Yellow	Flex 200, 250 Univ, 300, 200 Ext
PMP-020-045U	4642090	100 to 1000	1	Blue	Flex 1000, 1000, 1000 Ext
PMP-020-050E	4642100	0.5mL to 5mL	10	Green	5 mL
PMP-020-055R	4642110	1mL to 10mL	20	Red	10mL

Accessories

Catalogue No	Alt. No	Description
PMP-010-500G	9420400	Pipettor F-stand, 6 position