

Pipettors, electronic, single and multichannel,
Thermo Scientific Finpipette® Novus



04



12 channel

Catalogue No	Alt. No	Volume range, µL	Increments, µL	Colour code	Finntip
PMP-448-020L	46300100	1 to 10	0.1	Pink	10, 20, 50
PMP-448-040F	46300300	5 to 50	0.1	Yellow	250
PMP-448-060W	46300500	30 to 300	1.0	Orange	300, Flex 300

16 channel

Catalogue No	Alt. No	Volume range, µL	Increments, µL	Colour code	Finntip
PMP-448-070T	46300700	5 to 50	0.5	Turquoise	50



Pipettor tips, Thermo Scientific Finntip®



04



Accuracy and precision depend on the quality of the tips used for pipetting. Designed to optimise the performance of Thermo Scientific Finpipette®. Advanced moulding techniques and strict quality control ensure a perfect fit. Specialist tips, including the Finntip® 250 and 1,000 range, have a bevelled orifice which reduces shear force and thereby give added security to pipetting fragile items. Finntip® 10/20/50 are autoclavable.

Finntip® 10

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-107-090C	9400 310	0.2 to 10	Bag	-	1,000
PMP-107-110W	9400 300	0.2 to 10	Rack	96	960
PMP-107-113Q	9400 303	0.2 to 10	Sterile rack	96	960
PMP-107-121R	9400 326	0.2 to 10	Rack refill starter kit	192	2,688
PMP-107-122P	9400 37	0.2 to 10	Rack refill kit	192	3,840

Finntip Flex 10

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-107-500J	94060100	0.2 to 10	Rack	96	960
PMP-107-505W	94060103	0.2 to 10	Sterile rack	96	960
PMP-107-510G	94060116	0.2 to 10	Refill starter kit	96	2,016
PMP-107-515T	94060117	0.2 to 10	Refill kit	96	1,920
PMP-107-520D	94060120	0.2 to 10	Bag	-	1,000

Finntip® 20

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-108-010U	9400 610	0.2 to 20	Rack	384	3,840
PMP-108-020R	9400 613	0.2 to 20	Sterile rack	384	3,840
PMP-108-030X	9400 620	0.2 to 20	Bag	-	1,000

Finntip® 50

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-160-509X	9400 370	0.2 to 50	Rack	384	3,840
PMP-160-507S	9400 373	0.2 to 50	Sterile rack	384	3,840
PMP-160-505W	9400360	0.2 to 50	Bag	-	1,000

Finntip® Geno DNA, 10 to 250µL

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-103-100D	9405020	10 to 250	Bag	-	1,000
PMP-103-104S	9405030	10 to 250	Bulk	-	22,000
PMP-103-110A	9405120	10 to 250	Rack	96	960
PMP-103-114P	9405123	10 to 250	Sterile rack	96	960

Finntip® 200, extended 78mm

Catalogue No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
PMP-105-030J	9400100	5 to 200	Box	-	400
PMP-105-035W	9400130	5 to 200	Rack	96	960
PMP-105-040G	9400133	5 to 200	Sterile rack	96	960
PMP-105-045T	9400115	5 to 200	Sterile individually wrapped	-	250



entry continued