



Pipettor tips, epTIPS standard

96

eppendorf



High quality pipettor tips, packed loose in re-sealable ZipLock bags to protect from contaminants. PMP-890-521Q is the long 10mL top. (243mm)

Catalogue No	Alt. No	Volume range, μL	Colour code	Inner pack qty	Pack qty
PMP-890-501W	0030000811	0.1 to 10	Anthracite	500	1,000
PMP-890-503S	0030000838	0.1 to 20	Dark grey	500	1,000
PMP-890-505X	0030000854	0.5 to 20	Light grey	500	1,000
PMP-890-507K	0030000870	2 to 200	Yellow	500	1,000
PMP-890-509G	0030000897	20 to 300	Orange	500	1,000
PMP-890-511T	0030000919	50 to 1,000	Blue	500	1,000
PMP-890-513P	0030000935	50 to 1,250	Green	500	1,000
PMP-890-515L	0030000951	500 to 2,500	Red	100	500
PMP-890-517H	0030000978	100 to 5,000	Violet	100	500
PMP-890-519D	0030000765	1,000 to 10,000	Turquoise	100	200
PMP-890-521Q	0030000781	1,000 to 10,000	Turquoise	100	200



Pipettor tips, filter, epTIPS Dualfilter

96

eppendorf



ep Dualfilter T.I.P.S. consist of a unique two phase filter protection barrier made from hydrophobic polyethylene. The upper blue layer provides very effective protection against small aerosols and biomolecules (DNA). The lower layer acts as pre-filter and efficiently protects against larger aerosols, drops, splashes and larger amounts of liquids. The combination of both layers leads to a unique efficiency against aerosols contaminants. All ep Dualfilter T.I.P.S. are PCR* clean (free of human DNA, DNase, RNase, PCR* inhibitors) and sterile (to USP, DAB, Ph.Eur. standards) to satisfy the particular requirements of the molecular biology laboratory.

Catalogue No	Alt. No	Volume range, μL	Colour code	Pack qty
PMP-906-500U	0030077504	0.1 to 10	Anthracite	960
PMP-906-505K	0030077512	0.1 to 10	Dark grey	960
PMP-906-510R	0030077520	0.5 to 10	Light grey	960
PMP-906-515H	0030077539	2 to 20	Yellow	960
PMP-906-520X	0030077547	2 to 100	Yellow	960
PMP-906-525E	0030077555	2 to 200	Yellow	960
PMP-906-530L	0030077563	20 to 300	Orange	960
PMP-906-535B	0030077571	50 to 1,000	Blue	960
PMP-906-541G	0030077580	100 to 5,000	Violet	120
PMP-906-545V	0030077598	1,000 to 10,000	Turquoise	100

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche



Pipettor tips, epTIPS Biopur®

96

eppendorf



Biopur® tips are guaranteed sterile and free of pyrogens, RNase, DNA and ATP. They satisfy the most stringent demands in medicine, the pharmaceutical and food industries, PCR* technology, molecular biology and cytotechnology. epTIPS racks come in ready to use boxes and epTIPS singles come in individual blister strips.

epTIPS racks

Catalogue No	Alt. No	Volume range, μL	Colour code	Inner pack qty	Pack qty
PMP-900-502G	0030075005	0.1 to 20	Dark grey	96	480
PMP-900-503E	0030075021	2 to 200	Yellow	96	480
PMP-900-507T	0030075048	20 to 300	Orange	96	480
PMP-900-509P	0030075064	50 to 1,000	Blue	96	480
PMP-900-513B	0030075080	50 to 1,250	Green	96	480
PMP-900-515U	0030075102	500 to 2,500	Red	48	240
PMP-900-517Q	0030075145	1,000 to 10,000	Turquoise	24	120

epTIPS singles

Catalogue No	Alt. No	Volume range, μL	Colour code	Inner pack qty	Pack qty
PMP-900-505A	0030010035	2 to 200	Yellow	1	100
PMP-900-511F	0030010051	50 to 1,000	Blue	1	100

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche