

Pumps, vacuum and pressure, oil free

04


2^{yr}

Oil free pumps for degassing, filtration, desiccation, aspiration and tissue culture.

- Compact
- Economical
- Aqueous and chemical PTFE versions
- High flow versions for viscous and large well systems
- Tissue culture versions accessorised ready for use

Technical Specification - Specific

Aqueous version	PXW-800-040T	PXW-800-060N	PXW-800-070K	PXW-800-090E
Pressure, mbar	133	93	80	80
Flow rate, L/ min	18	28	38	83
Dimensions [l x w x h], mm	194 x 114 x 191	206 x 224 x 254	206 x 224 x 254	381 x 254 x 254
Electrical supply	230V, 50Hz	230V, 50Hz	230V, 50Hz	230V, 50Hz
Features	Vacuum regulator Vacuum gauge Inlet catchpot Pressure regulator Pressure gauge Exhaust catchpot	Vacuum regulator Vacuum gauge Inlet catchpot Pressure regulator Pressure gauge Exhaust catchpot	Vacuum regulator Vacuum gauge Inlet catchpot Pressure regulator Pressure gauge Exhaust catchpot	Vacuum regulator Vacuum gauge Inlet catchpot N/A N/A N/A

Chemical PTFE version	PXW-805-010V	PXW-805-020S
Pressure, mbar	133	106
Flow rate, L/ min	18	42
Dimensions [l x w x h], mm	254 x 254 x 175	254 x 254 x 175
Electrical supply	230V 50/60Hz	230V 50/60Hz
Features	Low noise Self-draining heads Gas ballast	Low noise Self-draining heads Gas ballast

Aqueous version

Catalogue No	Alt. No	
PXW-800-040T	2522	
PXW-800-060N	2534	
PXW-800-070K	2546	
PXW-800-090E	2567	

Chemical PTFE version

Catalogue No	Alt. No	
PXW-805-010V	2012	
PXW-805-020S	2015	

Pumps, vacuum and pressure, oil free, general purpose

50


2^{yr}

For general purpose applications including vacuum ovens, glove boxes, transfer chambers, rough drying and desiccation.

- Fast oil free pumping
- Able to replace oil pumps in many applications
- Smaller and lighter than oil pumps
- Economical
- Fast pumping at low pressure
- Vacuum regulator
- Vacuum gauge
- Inlet catchpot

Technical Specification - Specific

	PXW-810-010V	PXW-810-020S	PXW-810-030P	PXW-810-040M
Pressure, mbar	5	10	10	80
Speed range, L/ min	50	45	80	168
Dimensions [l x w x h], mm	439 x 165 x 267	439 x 165 x 267	432 x 191 x 305	432 x 191 x 305
Electrical supply	230V 50/60Hz	230V 50/60Hz	230V 50/60Hz	230V 50/60Hz

Catalogue No	Alt. No	
PXW-810-010V	2561	
PXW-810-020S	2562	
PXW-810-040M	2585	

