



Vacuum gauge, dial, Büchi

61



NEW

- Dial gauge, 0mbar to 1,000mbar
- Needle valve control of vacuum setting
- Retort stand mounting

Manual vacuum regulation gives cost-effective solution for control of vacuum. Woulff bottle accessory highly recommended to act as a buffer vessel and condensate trap.

Technical Specification

Accuracy	±20mbar
Control	1,000mbar to 0mbar
Net Weight, kg	4.5

Catalogue No	Alt. No	
EVF-950-532D	47291	

Accessory

Catalogue No	Alt. No	Description
VAC-850-500M	047170	Woulff bottle 850 with holder



Pumps, vacuum and pressure, diaphragm, oil free

03



For vacuum filtration, degassing, impregnation and gas transfer. Quiet and clean, continuous oil free pumping. All media contacting parts made from aluminium and Viton® (FPM). IP54 anti-splash protection rated and with thermal overload protection. All pumps fitted with handle, on-off switch, mains cable and rubber feet. Double diaphragm for high reliability and long life. Easy to service.

Technical Specification - Specific

	PYB-525-010A	PYB-525-040X
Flow rate, m³/hr	1.9	6.5
Vacuum, mbar	<80	<80
Compression, bar	2	2
Dimensions [w x d x h], mm	173 x 275 x 175	235 x 315 x 200
Mass, kg	6.3	15.5
Electrical supply	230 to 240V 50/60Hz	230V 50/60Hz

Catalogue No	Alt. No	
PYB-525-010A	ME2	
PYB-525-040X	ME8	

Spare parts

Catalogue No	Alt. No	Description
PYB-527-010L	696810	ME2 kit consisting of 1 x diaphragm and 2 x valves
PYB-527-020Y	696811	MZ2/ME4R kit consisting of 2 x diaphragms and 4 x valves
PYB-527-030F	696812	ME8 service kit comprising 4 x diaphragms, 2 x valves