

# Safety hand protection disposable

Gloves, nitrile, disposable, blue, N-DEX®, Class I Medical

S3



S



Manufactured from 100% nitrile and containing absolutely no natural rubber proteins. Available lightly powdered or powder free and in a variety of lengths and thicknesses, these gloves provide three times the puncture resistance of natural rubber or vinyl gloves of comparable thickness. The powder free versions feature a non-allergenic polymer gel coating on the inside that helps eliminate allergic reactions to gloves. All these gloves feature a rolled cuff. Chemical resistance data is not available for these products as they are Class 1 Medical devices to EN 455-1.

N-DEX® Medical Exam, powder free, thickness 0.10mm, length 240mm

Catalogue No	Size	Pack qty
SAR-587-010X	Extra small	100
SAR-587-020L	Small	100
SAR-587-030Y	Medium	100
SAR-587-040F	Large	100
SAR-587-050C	Extra large	100

N-DEX® Medical, powder free, thickness 0.15mm, long cuff, length 280mm

Catalogue No	Size	Pack qty
SAR-589-020W	Small	50
SAR-589-030T	Medium	50
SAR-589-040Q	Large	50
SAR-589-050N	Extra large	50

N-DEX® sterile, powder free, thickness 0.15mm, long cuff, length 280mm

Catalogue No	Size	Pack type	Pack qty
SAR-590-020A	Small	Pair	50
SAR-590-030U	Medium	Pair	50
SAR-590-040R	Large	Pair	50
SAR-590-050X	Extra large	Pair	50

Glove box holders, acrylic

41

Fisherbrand®



Wall mounted, clear acrylic holders for easy dispensing of gloves. Each holder includes a polyester foam insert to securely hold a range of glove box sizes. Fit ensures easy removal of gloves. Available in three sizes. Wall mount screws included.

Catalogue No	Description	Dimensions [w x d x h], mm
FB70943	1-box holder	141 x 260 x 110
FB70944	2-box holder	295 x 260 x 110
FB70945	3-box holder	295 x 385 x 110

Glove box holders, wire

41

Fisherbrand®



Epoxy coated steel wire racks which hold most major brands of glove boxes. Available as a single box holder or the three box holder. Wall mount screws included.

Catalogue No	Description	Dimensions [w x d x h], mm
FB70950	1-box holder	140 x 108 x 208
FB70951	3-box holder	250 x 108 x 454

Gloves, nylon, anti-static

04



Hard wearing, reversible glove. Made from anti-static nylon in filament form to give a comfortable, dust free glove which reduces the risk of sparks and hence explosions and fires in areas associated with inflammable chemicals. Washable without loss of anti-static properties.

Catalogue No	Size	Pack type	Pack qty
STP-220-010C	Small	Pair	1
STP-220-020W	Large	Pair	1