

## Scoops, stainless steel

41

**Fisherbrand®**

Cylindrical, without handle

Catalogue No	Length, mm	Width, mm	Pack qty
FB65003	70	30	5
FB65005	90	30	5
FB65007	105	30	5
FB65009	120	30	5

Cylindrical, with handle

Catalogue No	Length, mm	Width, mm
FB65011	190	50
FB65013	260	75
FB65015	320	100
FB65017	350	115



## Scoops, sampling, white polystyrene, non-sterile

03

**scienceware®**

These measures are ideal for retrieving small, specific volume samples of powders or liquids. Their unique flat bottom shape allows them to be used as weighing boats for additional sampling accuracy. Use a separate spoon for each material sampled to avoid cross contamination. All sizes are moulded from white, high impact polystyrene for durability and stiffness when sampling packed powders.

Catalogue No	Capacity, mL	Pack qty
DIS-910-010R	1.25	12
DIS-910-030L	2.5	12
DIS-910-050F	5.0	12
DIS-910-070W	8.0	12
DIS-910-090Q	20.0	12



## Scoops, sampling, white polystyrene, non-sterile

03

**scienceware®**

These sampling scoops are accurate measures for smaller quantities of dry or liquid materials. Accurate capacity is listed on each scoop and levelling the sample to the flat top edge ensures a repeatable and consistent volume measure. Use a separate spoon for each material sampled to avoid cross contamination. The flat bottom shape and lightweight allows for use as a weighing boat. These one-piece scoops are moulded from durable, white, high impact polystyrene.

Catalogue No	Capacity, mL	Pack qty
DIS-915-010K	63	12
DIS-915-030E	83	12
DIS-915-050V	125	12
DIS-915-070P	250	12



## Scoops, sampling, Sterileware™

03

**scienceware®****S**

Ideal for use with wide mouth containers and for samples up to 250mL. White, high impact polystyrene, individually wrapped.

Catalogue No	Alt. No	Capacity, mL	Dimensions [w x d x h], mm	Length, mm	Pack qty
DIS-330-030W	H36902-0000	60	25 x 76 x 38	178	100
DIS-330-050Q	H36904-0000	125	25 x 89 x 44	190	100
DIS-330-070K	H36906-0000	250	51 x 127 x 51	216	100



## Spoons, sampling, Sterileware™

03

**scienceware®****S**

The spoons come in four sizes for precise measurement and are ideal for small amounts of sample material. The spoons have a long handle for easy reaching and to minimise contact with sample material. Individually wrapped.

Catalogue No	Alt. No	Capacity, mL	Length, mm	Inner pack qty	Pack qty
DIS-345-030Y	H36940-0000	1.25	178	1	200
DIS-345-050C	H36942-0000	2.5	178	1	200
DIS-345-070T	H36944-0000	5.0	178	1	200
DIS-345-090N	H36946-0000	15.0	178	1	200

