

Fisherbrand™ Cell Spreaders

L-Shaped Cell Spreaders

Ideal for even spreading of bacterial cultures


- Disposable
- Smooth spreading surface with slight upward angle reduces the risk of damaging agar interface
- L-shape allows application over the entire Petri dish or plate surface
- Presterilised blue polypropylene eliminates the need for flaming or autoclaving

Specifications

Length	145.1mm
Colour	Blue
Sterile	Yes
For Use With (Application)	Ideal for even spreading of bacterial cultures, smooth spreading surface with slight upward angle reduces the risk of damaging agar interface
Shape	L-shaped
Material	Polypropylene
Disposable	Yes
For Use With (Equipment)	Agar culture plates

Ordering Information

Cat. No.	Description	Pack qty
15615467	L-shaped cell spreader, subpack of 10	500
15625467	L-shaped cell spreader, individually wrapped	100



Triangular Cell Spreaders

Nonbreakable and reusable polypropylene spreaders available in two sizes: 30 and 60mm bar widths

Specifications

Autoclavable	Yes
Material	Polypropylene
Sterile	Yes
Shape	Triangle
For Use With (Equipment)	Use 30mm wide spreaders for prepared biplates; 60mm for 100 to 150mm prepared plates

Ordering Information

Cat. No.	Description	Pack qty
15635467	Triangular cell spreader blue 30mm, individually wrapped	25
15645467	Triangular cell spreader yellow 60mm, individually wrapped	25

