

Spatula weighing set, micro

03

**a**

Set of 6 micro size double-ended tools and a pair of tweezers for accurate handling of minute quantities. Each tool has two different end shapes, doubling the selection of implements to better match the users' needs. Tool ends include two small spoons for transferring and collecting both dry and liquid materials. There are also eight flat spatulas of various widths, lengths and end shapes, that are ideal for mixing, scraping and spreading. The tweezers have a serrated inside surface for secure sample handling and grooved handles for better gripping. Approximate length of instruments 140mm.

All tools are manufactured from polished stainless steel and can be steam autoclaved at 121°C. When not in use the tools can be stored inside the carrying case provided.

Catalogue No	Alt. No
SMJ-200-010P	F36705-0000




● Spatulas, palette knives

03

With parallel sided, flexible stainless steel blade set in plastic handle. Sterilisable.

Catalogue No	Blade length, mm	Blade width, mm
SMK-401-030A	102	20
SMK-401-050R	152	25
SMK-401-070L	203	30
SMK-401-090F	250	40
SMK-401-110C	300	45



● Spatulas, palette knives

03

With tapered, flexible stainless steel blade set in plastic handle. Sterilisable.

Catalogue No	Blade length, mm	Blade width, mm
SMK-461-030E	75	15
SMK-461-050V	100	20
SMK-461-070P	150	25

Spatula, filling knife

03

With stainless steel blade set in plastic handle. Sterilisable.

Catalogue No	Blade length, mm	Blade width, mm
SMK-510-030Y	130	50



Spatula, micro

46

Fisherbrand®

Colour coded for easy identification these flat ended micro spatulas have stainless steel blades and polystyrene handles. Overall length 60mm and blade widths 0.1mm, 0.3mm, 0.5mm and 1mm, making them ideal for handling small samples with ease.

Catalogue No	Pack qty
FB70365	4

