

Fisherbrand Comfort Nitrile gloves



Comfort Nitrile's sizes

	Size XS (5-6)	Size S (6.5 - 7)	Size M (7.5 - 8)	Size L (8.5 - 9)	Size XL (9.5 - 10)
Product Style / ID	15622367	15632367	15642367	15652367	15662367

Attributes

Product Category /Structure:	Disposable Glove
Material:	Acrylonitrile-butadiene (Nitrile)
Cuff:	Beaded
Powder Mass /Residue:	Max. 2.0 mg/glove
Texture:	Finger Tip
Color:	Periwinkle Blue
Special Attributes /Certifications:	NA
Product Class /Grade:	Industrial

Glove Configuration /Design:	Ambidextrous
Expiry date /Shelf Life:	3 year
Protein Level:	NA, not made from NRL
Product Style:	General Disposable
Lining:	NA
Product conformance:	PPE, Category III risk EN 420 EN ISO 374/ Type B Certified for food contact

Length

Typical Length (mm / in)	240 / 9.45
--------------------------	-------------------

Minimum Length (mm / in)	Min 230 / 9.06
--------------------------	-----------------------

Dimensions

	Size XS (5-6)	Size S (6.5 - 7)	Size M (7.5 - 8)	Size L (8.5 - 9)	Size XL (9.5 - 10)
Palm Width (mm)	70 ± 10	80 ± 10	95 ± 10	110 ± 10	120 ± 10

Thickness

Location of Thickness Measurements	Typical Single Wall	Minimal Single wall
Cuff (at 25 mm from cuff end)	0.06 mm / 2.4 mil	0.04 mm/ 1.6 mil
Palm (at center of palm)	0.07 mm / 2.8 mil	0.05 mm/ 2.0 mil
Finger (at 13.3mm from extreme tip)	0.09 mm/ 3.5 mil	0.07 mm / 2.8 mil

Physical Properties

	Typical Before Aging	Typical After Aging	Minimum Before Aging	Minimum After Aging
Tensile Strength (Mpa)	20	18	Min 14	Min 14
Elongation (%)	500	400	Min 500	Min 400

Criteria

	AQL	Inspection Level ASTM
Freedom from Holes	1.5	G-I
Dimensions	4.0	S-2
Physical Properties	4.0	S-2
Visual Inspections (Major)	2.5	N/A
Visual Inspections (Minor)	4.0	N/A
Powder Free Residue	N/A	N = 5

Fisherbrand Comfort Nitrile

	Inches (L x W x H)	mm (L x W x H)
Inner Carton Dimensions	9.64 x 5.0 x 3.15	245 x 127 x 80
Outer Case Dimensions	16.54 x 10.24 x 10.04	420 x 260 x 255


Cube (ft)	Dispenser Box per Case	Gloves per Dispenser Box by Weight
0.98	10	XS-L – 200 / XL - 180

Product Style / ID (SKU)	Size	Inner Carton UPC Code	Outer Case UPC Code
15622367	Size XS (5-6)	15622367	NA
15632367	Size S (6.5 - 7)	15632367	NA
15642367	Size M (7.5 - 8)	15642367	NA
15652367	Size L (8.5 - 9)	15652367	NA
15662367	Size XL (9.5 - 10)	15662367	NA

CAUTION: THIS PRODUCT TO BE STORED IN CARTON DURING TRANSPORTATION

Permeation Performance Levels Per EN ISO 374:2016

Level	Breakthrough Time (min)	Chemical	Levels
1	>10		
2	>30		
3	>60	n-Heptane	6
4	>120	40% Sodium Hydroxide	6
5	>240	Formaldehyde	6
6	>480		



© 2017 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500