


Spatulas, polypropylene, tapered

02

Azlon

Catalogue No	Alt. No	Length, mm	Width, mm	Pack qty
SMJ-790-011C	SWY060	100	17	10
SMJ-790-031T	SWY062	150	26	10
SMJ-790-051N	SWY064	200	28	10
SMJ-790-071H	SWY066	250	30	10



Spatula, sampling, Sterileware™

03

Ideal for every type of powdered or granulated material. This 229mm long spatula accommodates twice the volume of a standard reusable spatula, facilitating sampling from large containers, bags or straight from production lines. The sturdy edges may be used to break up any caked material while sampling. Individually wrapped.

Catalogue No	Alt. No	Length, mm	Pack qty
DIS-340-030P	H36930-0000	229	100



See also Scoops pages 1,252 to 1,254

Spatula, wooden applicator

03

Catalogue No	Pack qty
DIS-295-037N	864

