

Stirrers, magnetic, unheated, mini, Topolino / Topolino mobil



60

NEW



Topolino:

- Durable, brushless motor
- Continuously adjustable speed range
- High magnetic connection

Topolino mobil:

Same features as the Topolino, plus:

- Portable operation with long operating time (8 to 12hr)
- Short charging time (2 to 3hr)
- Standard replaceable AA rechargeable batteries
- Optional power mode:
 - Mains-free with standard batteries
 - With supplied mains adapter (without batteries)
 - Combined mains/battery operation (with batteries fitted)



Technical Specification - Specific

| | SWT-840-010F | SWT-840-020C |
|----------------------------|---------------|------------------------------|
| Dimensions [w x d x h], mm | 95 x 110 x 40 | 140 diameter x 40 |
| Weight, g | 300 | 320 |
| Electrical supply | 230V 50/60Hz | 230V 50/60Hz AA batteries |

Technical Specification - General

| | |
|---------------------------------------|------------------------------|
| Stirring quantity H ₂ O, L | 0.25 (max) |
| Motor rating, input output, W | 1.1 / 0.8 |
| Speed range, rpm | 300 to 1.800 |
| Max. magnetic bar [L x Ø], mm | 40 x 6 |
| Material | Polypropylene set up plate |
| Plate diameter, mm | 80 |
| Ambient temperature, °C | 5 to 40 (permissible) |
| Relative humidity | 80% (permissible) |
| Protection class | IP 21 (acc. to DIN EN 60529) |

| Catalogue No | Alt. No | Model No |
|--------------|---------|----------------|
| SWT-840-010F | 3368000 | Topolino |
| SWT-840-020C | 3381300 | Topolino mobil |

Stirrer, magnetic, unheated, air operated



Manufactured from high impact polystyrene, this very effective air powered magnetic stirrer is simple to use and very economical to run. The only additional requirement is a low pressure air supply to power the magnetic stirrer, making it ideal for environments where sparking hazards from electrical sources must be eliminated. This compact stirrer operates using only 3psi to 8psi (21kPa to 55kPa) air pressure and is ideal for stirring vessels up to 1L. With tubing to connect to an air supply.



| Catalogue No | Alt. No | Height, mm | Diameter, mm |
|--------------|-------------|------------|--------------|
| SWT-165-010R | F37002-0000 | 20 | 100 |

Stirrers, magnetic, unheated, lab disc IKAMAG®



60

- Ultra flat
- Advanced magnetic coil with no moving parts
- Wear-free
- Reverses direction of rotation every 30s to ensure better mixing



Reliable, strong mixers with great designs that bring colour to the laboratory.

Technical Specification - General

| | |
|---------------------------------------|------------------|
| Stirring quantity H ₂ O, L | 0.8 (max) |
| Motor rating, input output, W | 5 / 3 |
| Speed range, rpm | 15 to 1,500 |
| Magnetic stirrer bar, mm | 30 x 8 (maximum) |
| Plate diameter, mm | 90 |
| Dimensions [w x d x h], mm | 114 x 161 x 12 |
| Mass, kg | 0.3 |
| Electrical supply | 230V 50Hz |

| Catalogue No | Alt. No | Description |
|--------------|---------|------------------------|
| SWT-865-010H | 3362300 | Lab disc Sunny Side Up |
| SWT-865-030B | 3277700 | Lab disc White |
| SWT-865-150X | 3578000 | Lab disc Orange |
| SWT-865-155E | 3578600 | Lab disc Brown |
| SWT-865-160L | 3578200 | Lab disc Green |
| SWT-865-165B | 3578400 | Lab disc Pink |