

Stirrers magnetic unheated

Stirrer, magnetic, unheated, KMO 2 basic, IKAMAG®

60

2^{yr}



- Small, powerful, unheated magnetic stirrer
- Infinitely variable speed from 0 to 1,100rpm
- Strong magnetic field
- Motor with optoelectronic speed control
- Stainless steel casing facilitates cleaning and sterilisation
- Includes M 10 thread for H16 V support rod

Technical Specification

Stirring quantity H ₂ O, L	5
Motor rating, input output, W	14 / 4
Speed range, rpm	0 to 1,100
Plate dimensions, mm	140 x 120
Dimensions [w x d x h], mm	140 x 200 x 75
Mass, kg	1.4
Electrical supply	230V 50/60Hz
Protection class	IP21

Catalogue No	
SWT-385-010K	

Accessory

Catalogue No	Alt. No	Description
SWT-565-510L	H 16 V	H 16 V support rod

IKA®

Stirrers, magnetic, unheated, digital, REO basic IKAMAG®

60

IKA®

NEW



- Digital speed display
- Non-locking, electronically controlled motor
- Constant speed even during changes in load
- Infinitely variable speed
- Incl. protection cover H 101

Technical Specification

Stirring quantity H ₂ O, L	20
Speed display	digital
Speed range, rpm	0 to 1,200
Max. magnetic bar [L x Ø], mm	80 x 10
Material	stainless steel set up plate
Plate diameter, mm	Ø 135
Dimensions [w x d x h], mm	160 x 270 x 85
Weight, kg	2.4
Ambient temperature, °C	5 to 40 (permissible)
Relative humidity	80% (permissible)
Protection class	IP 42 (acc. to DIN EN 60529)
Electrical supply	230V 50/60Hz

Catalogue No	Alt. No	
SWT-388-010P	3384200	

Accessories

Catalogue No	Alt. No	Description
SWT-755-500F	1572000	Stir bar 10mm x 6mm, round PTFE coated
SWT-755-505S	1572100	Stir bar 15mm x 6mm, round PTFE coated
SWT-755-510C	1572200	Stir bar 20mm x 8mm, round PTFE coated
SWT-755-515P	1572300	Stir bar 25mm x 8mm, round PTFE coated
SWT-755-520W	1572400	Stir bar 30mm x 8mm, round PTFE coated
SWT-755-525M	1572500	Stir bar 40mm x 8mm, round PTFE coated
SWT-755-530T	1572600	Stir bar 50mm x 8mm, round PTFE coated
SWT-755-535J	1572800	Stir bar 80mm x 8mm, round PTFE coated
SWT-755-540Q	0793300	Stir bar 110mm x 27mm, round PTFE coated
SWT-755-545G	1129000	Stir bar 155mm x 27mm, round PTFE coated
SWT-755-550N	0356600	Stir bar, triangular, 25mm length, PTFE coated
SWT-755-560K	0356500	Stir bar, triangular, 40mm length, PTFE coated