

## Stirrer hotplate, magnetic, RH basic 2 IKAMAG® 04

2<sup>yr</sup>

- Economic magnetic stirrer with stainless steel heating plate
- Fixed safety circuit 400°C
- Soft-start stirring motor

**IKA®**



### Technical Specification

Stirring quantity H <sub>2</sub> O, L	10
Power, W	15 (input) / 2 (output)
Display	Scale (0 to 6)
Speed range, rpm	100 to 2,000
Output	400W
Heating rate, K/min	3
Temperature range, °C	Ambient to 320°C
Diameter, mm	125 (hotplate)
Dimensions [w x d x h], mm	168 x 220 x 105
Mass, kg	2.4
Ambient temperature range	5 to 40
Humidity, %RH	80
Protection class	IP 21
Electrical supply	230V, 50/60Hz

Catalogue No	
SWT-145-010Y	

### Accessories and spares

Catalogue No	Alt. No	Description
TEH-270-010Q	ETS-D5	Contact thermometer, with pH measurement facility
TEH-270-020N	ETS-D6	Contact thermometer
SWT-145-500X	H15	1.5L bath attachment, 140mm inner dia. x 125mm, for tempering 0.5 and 1L flasks
SWT-145-505E	H28	1L bath attachment, 140mm inner dia. x 70mm, suitable for sand bath. (max. temperature 350°C)
SWT-755-500F	1572000	Stir bar 10mm x 6mm, round PTFE coated
SWT-755-505S	1572100	Stir bar 15mm x 6mm, round PTFE coated
SWT-755-510C	1572200	Stir bar 20mm x 8mm, round PTFE coated
SWT-755-515P	1572300	Stir bar 25mm x 8mm, round PTFE coated
SWT-755-520W	1572400	Stir bar 30mm x 8mm, round PTFE coated
SWT-755-525M	1572500	Stir bar 40mm x 8mm, round PTFE coated
SWT-755-530T	1572600	Stir bar 50mm x 8mm, round PTFE coated
SWT-755-535J	1572800	Stir bar 80mm x 8mm, round PTFE coated
SWT-755-540Q	0793300	Stir bar 110mm x 27mm, round PTFE coated
SWT-755-545G	1129000	Stir bar 155mm x 27mm, round PTFE coated
SWT-755-550N	0356600	Stir bar, triangular, 25mm length, PTFE coated
SWT-755-560K	0356500	Stir bar, triangular, 40mm length, PTFE coated
SWT-755-565A	1358600	Set of stir bars, length 10mm to 80mm
SWT-765-570A	1293100	Stir bar remover