

# Stirrers magnetic followers

Magnetic followers, PTFE, single sided star head

**Fisherbrand®**

41



Catalogue No	Length, mm	Diameter, mm	Pack qty
FB55660	8	10	5
FB55661	12	14	5
FB55662	13	17	5
FB55663	15	25	5
FB55664	17	40	1
FB55665	17	60	1

Magnetic followers, PTFE, micro

**Fisherbrand®**

41



For stirring in test tubes, small size beakers or other containers.

Supplied in re-useable snap cap vials.


Catalogue No	Length, mm	Diameter, mm	Pack qty
FB55605	2	2	10
FB55606	5	2	10
FB55607	8	1.5	10
FB55609	7	2	10
FB55611	10	3	10

Magnetic followers, PTFE, Saturn® Spinbar

**Fisherbrand®**

41

**a**



Manufactured from inert and chemically resistant PTFE the Saturn® Spinbar is autoclavable at 121°C. The Saturn® has a large bulb at the centre which lifts the stir bar at the ends whilst spinning, and thereby reduces the surface contact area allowing the bar to spin with decreased drag. Available in 2 lengths to match a wide variety of vessels, with both round and flat bottoms.

Catalogue No	Length, mm
FB55620	40
FB55621	50

Magnetic followers, PTFE, cross

**Fisherbrand®**

41

**NEW a**



A stable, general purpose stirrer. Provides a deep, quiet vortex even at relatively high speeds.

Catalogue No	Length, mm	Diameter, mm	Pack qty
FB58190	10	5	5
FB58191	20	8	5
FB58192	25	9	5
FB58193	30	10	5
FB58194	38	11	5

Magnetic followers, PTFE, coloured, double ended

**Fisherbrand®**

41

**a**



Ideal for colour coded experiments. The double paddle action gives efficient stirring plus high stability.

Magnetic followers, PTFE, oval

**Fisherbrand®**

41



Catalogue No	Length, mm	Diameter, mm	Pack qty
FB55630	10	5	5
FB55631	15	6	5
FB55632	20	10	5
FB55633	25	12	5
FB55634	30	16	5
FB55635	35	16	5
FB55636	40	20	1
FB55637	50	20	1
FB55638	64	20	1
FB55639	70	20	1

Catalogue No	Length, mm	Diameter, mm	Colour
FB58160	35	8	Red
FB58161	35	8	Blue
FB58162	35	8	Yellow
FB58163	35	8	White
FB58165	55	8	Red
FB58166	55	8	Blue
FB58167	55	8	Yellow
FB58168	55	8	White