

Stirrers overhead

Stirrers, overhead, small,
RW 11 basic 'Lab Egg'

60

IKA®

Stirrer, overhead, RW 14

04

IKA®

2^{yr}



- Glass housing resistant to chemicals
- Max. stirring quantity 2L (H₂O)
- Incl. paddle stirrer R 1001 and extension arm

Small-sized stirrer available in four colours: transparent, off white, creamy blue, and salmon pink.

Technical Specification - General

Stirring quantity H ₂ O, L	20
Speed range, rpm	0 to 2,000
Viscosity, mPas	100 (maximum)
Motor rating, input output, W	8 / 1
Mass, kg	0.39
Ambient temperature range	5 to 40°C
Electrical supply	230V 50/60Hz

Catalogue No	Colour
1 SWR-800-010H	Transparent
2 SWR-800-020E	Off white
3 SWR-800-050S	Creamy blue
SWR-800-060P	Salmon pink

Accessories and spares

Catalogue No	Description
SWR-850-505X	R 103 stand
SWR-850-510V	Propeller stirrer
SWR-850-520S	Vane type stirrer

NEW



- Constant speed
- Infinitely adjustable speed
- Easy to operate
- Slimline
- Non-locking, overload capabilities

Quiet, economical laboratory stirrer with infinitely adjustable speed.

For stirring substances of low to medium viscosity.

Technical Specification

Stirring quantity H ₂ O, L	8
Viscosity, mPas	10,000 (maximum)
Motor rating, input output, W	28.5 / 17
Output	17W (at stirring shaft)
Max. On time	100%
Torque, Ncm	8 (maximum at chuck)
Speed range, rpm	100 to 2,000
Speed display	scale
Chuck range, mm	0.5 to 8
Diameter, mm	13
Length, mm	160 (extension arm)
Dimensions (w x d x h), mm	70 x 176 x 197 (without extension arm)
Weight, kg	2
Ambient temperature, °C	5 to 40 (permissible)
Relative humidity	80% (permissible)
Protection class	IP 20 (acc. to DIN EN 60529)
Electrical supply	220-240V 50/60Hz

Catalogue No	Alt. No	
SWR-335-010T	3331400	

Accessories

Catalogue No	Alt. No	Description
SWR-535-515T	3160000	Plate stand, 200mm (w) x 316mm (d) x 560mm (h)
TWK-614-501Y	3160100	Plate stand, 200mm x 316mm x 800mm (w x d x h)
TWK-614-515U	3160200	Plate stand, 200mm x 316mm x 1,000mm (w x d x h)
TWK-614-505A	2657700	Bosshead
TWK-614-510H	3008600	Strap clamp
SWR-345-505Y	2603000	Shaft protection, Plexiglass, length adjustment 190mm/310mm
SWR-335-500C	2604000	Support holder R301.1
TWK-614-525R	2336000	Flexible coupling for attaching glass stirring rods
SWR-535-525Q	2022300	Anchor stirrer, 350mm (l) x 45mm (blade width)
SWR-350-560A	0741000	Propellor stirrer, 4-bladed, 350mm (l) x 50mm (blade width)
SWR-350-565N	0757600	Paddle stirrer R1373, 350mm (l) x 70mm (blade width)