



Continuation Stirrers, overhead EUROSTAR

60



2 yr

Stirring rod compatibility

	SWR-345-010M	SWR-345-020J	SWR-345-030G	SWR-345-040D	SWR-345-050A	SWR-345-060U
SWR-350-560A	X	X	X	X	X	X
SWR-350-570U			X	X	X	X
SWR-360-500L	X	X	X	X	X	X
SWR-360-505B	X	X	X	X	X	X
SWR-360-530C	X	X	X	X	X	X
SWR-360-510Y	X	X	X	X	X	X
SWR-345-500S	X	X	X	X	X	X
SWR-360-515V						
SWR-350-585H	X	X	X	X		
SWR-350-590X		X	X	X		
SWR-350-592K	X	X	X	X		
SWR-350-610L	X	X	X	X	X	X
SWR-350-595E	X	X	X	X	X	X
SWR-350-600X		X	X	X	X	X
SWR-350-605E					X	X
SWR-535-525Q	X	X	X	X	X	X
SWR-345-525C		X	X	X	X	X
SWR-535-530A	X	X	X	X	X	X
SWR-350-615B	X	X	X	X	X	X

Catalogue No	Description
SWR-345-010M	Digital
SWR-345-020J	Power basic
SWR-345-030G	Power control-visc
SWR-345-040D	Power control-visc, P1
SWR-345-050A	Power control-visc, P4
SWR-345-060U	Power control-visc, P7

Accessories for EUROSTAR digital

Catalogue No	Alt. No	Description
SWR-535-515T	3160000	Plate stand, 200mm x 316mm x 560mm (w x d x h)
TWK-614-501Y	3160100	Plate stand, 200mm x 316mm x 800mm (w x d x h)
TWK-614-515U	3160200	Plate stand, 200mm x 316mm x 1,000mm (w x d x h)
TWK-614-505A	2657700	Bosshead
TWK-614-510H	3008600	Strap clamp
SWR-345-505Y	2603000	Shaft protection, Plexiglass, length adjustment 190mm/310mm

Accessories for EUROSTAR power basic and control-visc

Catalogue No	Alt. No	Description
MPR-350-525N	1412000	H stand, rod height 1,010mm
SWR-535-500J	1412100	Telescopic stand, 620mm to 1,010mm
SWR-535-510G	2664000	Bosshead clamp
TWK-725-505C	3159000	Strap clamp
SWR-345-505Y	2603000	Shaft protection, Plexiglass, length adjustment 190mm/310mm
SWR-345-530J	2970000	labworldsoft® software (for control-visc stirrers only)

Extra accessory for EUROSTAR power control-visc

Catalogue No	Description
SGM-187-535R	PC 1.5m cable

Stirrer rods

Catalogue No	Alt. No	Description
SWR-350-560A	0741000	Propellor stirrer, 4-bladed, 350mm (l) x 50mm (blade width)
SWR-350-570U	0741300	Propellor stirrer, 4-bladed, 540mm (l) x 100mm (blade width)
SWR-360-500L	R1381	Propellor stirrer, 3-bladed, 350mm (l) x 45mm (blade width)
SWR-360-505B	R1382	Propellor stirrer, 3-bladed, 350mm (l) x 55mm (blade width)
SWR-360-530C	R1389	Propellor stirrer, PTFE coated, 3-bladed, 350mm (l) x 75mm (blade width)
SWR-360-510Y	R1311	Turbine stirrer, 350mm (l) x 30mm (blade width)
SWR-345-500S	R1312	Turbine stirrer, 350mm (l) x 50mm (blade width)
SWR-360-515V	2333100	Turbine stirrer, 400mm (l) x 70mm (blade width)
SWR-350-585H	0513500	Dissolver stirrer, 350mm (l) x 80mm (blade width)
SWR-350-590X	2387900	Dissolver stirrer, 350mm (l) x 100mm (blade width)
SWR-350-592K	R1303	Dissolver stirrer, 350mm (l) x 42mm (blade width)
SWR-350-610L	R1352	Centrifugal stirrer, 350mm (l) x 15/60mm (blade width)
SWR-350-600X	0757700	Paddle stirrer, 550mm (l) x 70mm (blade width)
SWR-350-605E	0757800	Paddle stirrer, 550mm (l) x 150mm (blade width)
SWR-535-525Q	2022300	Anchor stirrer, 350mm (l) x 45mm (blade width)
SWR-345-525C	2022400	Anchor stirrer, 350mm (l) x 90mm (blade width)
SWR-535-530A	R1332	Anchor stirrer, PTFE coated, 350mm (l) x 60mm (blade width)
SWR-350-615B	R1335	Kneading stirrer, 350mm (l) x 45mm (blade width)