

Stoppers, polyurethane foam, white

41

Fisherbrand®
a

Convenient method of 'stopping' sterile media flasks and vials. A reliable alternative to the cotton plugs normally used for this purpose. During storage of media, it is preferable to cover the stoppers (i.e. over the neck of the media container) with aluminium foil to avoid dust contamination on removal. Shrinkage will occur after continued use following autoclaving, but this is dependent purely on application. To eliminate the possibility of obstructing the flow of air into the media container, care should be taken not to crease the stoppers on fitting.

Catalogue No	Diameter, mm	Length, mm	Pack qty
FB68836	25	50	100
FB68837	38	50	100
FB68838	50	50	100



Stoppers, rubber, turnover

41

Fisherbrand®

Turnover stoppers create a very effective double seal for narrow aperture containers. The external serrations create an excellent internal seal, with the flexible turnover sleeve creating a secure seal on the container's outer neck. The head is easily pierced with a hypodermic needle.

Red

Catalogue No	Plug diameter, mm	Pack qty
FB57872	6.5	10
FB57873	8	10
FB57874	9.5	10
FB57875	11	10
FB57876	12.5	10
FB57877	14	10
FB57878	16	10
FB57879	17.5	10
FB57880	19	10
FB57881	20.5	10
FB57882	22	10
FB57883	24	10
FB57884	25.5	10



White

Catalogue No	Plug diameter, mm	Pack qty
FB68680	6.5	10
FB68681	8	10
FB68682	9.5	10
FB68683	11	10
FB68684	12.5	10
FB68685	14	10
FB68686	16	10
FB68687	17.5	10
FB68688	19	10
FB68689	20.5	10
FB68690	22	10
FB68691	24	10
FB68692	25.5	10

Blue

Catalogue No	Plug diameter, mm	Pack qty
FB68824	12.5	10
FB68827	17.5	10
FB68829	20.5	10

Stoppers, silicone, plug

41

Fisherbrand®
a

FDA approved silicone rubber to USP Class VI and BgVV, these pull stoppers are non-toxic and autoclavable. Externally serrated to create an effective internal seal with an elliptical flange to allow efficient removal of stopper from aperture. Easily pierced with a hypodermic needle they will fit standard tube sizes.

Catalogue No	Plug diameter, mm	Description	Pack qty
FB57895	10.7	Fits tube diameter 9mm	10
FB57896	14.6	Fits tube diameter 14mm	10
FB57897	16.5	Fits tube diameter 16mm	10
FB57898	23	Fits tube diameter 22mm	10

