

Thermometers liquid filled



Thermometers, white backed, red spirit filled,
total immersion

Fisherbrand®

41

General purpose/stirring thermometers with amber graduations, reinforced bulb for stirring. Non-roll fitting.

Catalogue No	Range, °C	Subdivisions, °C	Length, mm
FB58541	-10 to +50	0.5	305
FB58543	-20 to +110	1	305
FB58545	-20 to +150	1	305
FB58553	-10 to +50	1	155
FB58555	-10 to +110	1	155



Thermometers, white backed, red spirit filled,
76mm, immersion

Fisherbrand®

41

General purpose/stirring thermometer with amber graduations, reinforced bulb for stirring. Non-roll fitting.

Catalogue No	Range, °C	Subdivisions, °C	Length, mm
FB58547	-10 to +50	0.5	305
FB58549	-10 to +110	1	305
FB58551	-10 to +150	1	305



Thermometers, mercury free, blue spirit filled,
Lo-Tox™



33

Catalogue No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
THL-460-010D	44/812/8	-10 to 210	1	405
THL-460-030U	44/809/8	-35 to 50	1	305
THL-460-040R	44/808/8	-10 to 50	0.5	305
THL-460-050X	44/810/8	-20 to 110	1	305
THL-460-060L	44/811/11	-20 to 150	1	305
THL-460-110W	44/830/7	-10 to 110	1	115
THL-460-120T	44/850/5	-10 to 50	1	155

Catalogue No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
THL-460-020A	44/804/8	-10 to 260	1	405
THL-460-070Y	44/801/8	-35 to 50	1	305
THL-460-080F	44/800/8	-10 to 50	0.5	305
THL-460-090C	44/802/8	-20 to 110	1	305
THL-460-100C	44/803/8	-20 to 150	1	305



Thermometers, mercury free, blue spirit filled,
Lo-Tox™



33

General purpose/stirring thermometer with a low toxic, environmentally friendly blue filling, which is easy to read against the white background. Each thermometer is manufactured to traceable standards, and comes with a blue non-roll top and safety suspension ring.

Fractional degree
Total immersion

Catalogue No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
THL-465-010T	44/821/0	-1 to 51	0.1	405
THL-465-020Q	44/823/0	-1 to 51	0.2	405
THL-465-030N	44/824/0	-1 to 101	0.2	405

76mm immersion

Catalogue No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
THL-465-040K	44/820/0	-1 to 51	0.1	405
THL-465-050H	44/822/0	-1 to 51	0.2	405
THL-465-060E	44/824/0	-1 to 101	0.2	405