

# Timers

## Timer, 3 channel alarm, Traceable® 50 **Fisherbrand®**



- Adjustable alarm volume for quiet or noisy areas
- Program times from 1s to 100hr with memory
- Accuracy 0.001%
- Stopwatch
- Time-of-day clock
- 'Continuous alarm' setting ensures alarm is never missed
- Dimensions: 200mm x 2.5mm diameter, overall length: 256mm, dial: 43mm diameter
- With batteries
- Supplied with Traceable® certificate ISO 17025 A2LA to NIST

Catalogue No	
FB70228	

### Spare part

Catalogue No	Description	Pack qty
BMT-680-040W	Batteries alkaline, AAA (LR03)	4

## Timer, 4 channel alarm, Traceable® 50 **Fisherbrand®**



- Easy viewing digits are 19mm high
- Count down/up timing on four separate channels
- Program times from 1s to 100hr with memory
- Accuracy 0.01%
- Finger-sized keys are easy to use
- Displays time of day
- With stand and magnet
- Dimensions: 69mm x 38mm x 12.5mm
- With batteries
- Timing capacity 100h to 1s
- Supplied with Traceable® certificate ISO 17025 A2LA to NIST

Catalogue No	
FB70230	

### Spare part

Catalogue No	Description	Pack qty
BMD-140-011E	Battery alkaline 130mAh 1.5V LR44	4

## Timer, 3 line, Traceable® 50 **Fisherbrand®**

NEW



- Clock, count up/down, alarm, stopwatch
- Displays 3 different countdown/alarm times simultaneously
- Three unique channels may be started individually or concurrently
- Memories return the display to the previously programmed time
- Perfect unit to dedicate to a group of repetitive tests
- 6 digit LCD 9.5mm high display
- Countdown: 1s increments from 99hr, 59min, 59s
- Alarm: sounds at zero "Time's Up" is displayed, timer begins counting up
- Accuracy: 0.01%.
- Chemical-resistant ABS plastic case, magnetic back, flip-open easel for standing
- Supplied with 1½ year continuous-use silver-oxide battery
- Size: 51mm x 83mm x 13mm; weight: 85g
- Supplied with Traceable® certificate ISO 17025 A2LA to NIST

Catalogue No	
FB70295	

## Timer, high volume, Traceable® 50 **Fisherbrand®**

NEW



- Counts down and alarms
- 4 digit LCD 19mm high display shows the time remaining in minutes and seconds
- Countdown: 1s increments from 99m, 59s
- Automatic memory resets display to last programmed time, ideal for repetitive tests
- High volume alarm: sounds at zero for 30s
- Accuracy: 0.01%
- Magnetic back, flip-open easel for standing, a spring fastener for clipping to a lab coat
- Supplied with a 2 year, continuous-use, AAA battery
- Size: 63mm x 63mm x 13mm; weight: 57g
- Supplied with Traceable® certificate ISO 17025 A2LA to NIST

Catalogue No	
FB70325	