

## Test tubes, borosilicate glass, round bottom, screw thread with polypropylene cap, disposable 40

**Fisherbrand®**

Ideal for tissue culture, bacteriology, clinical chemistry, blood typing and cross matching procedures.

With natural, linerless, polypropylene screw cap and marking spot. Tubes packed in shrink wrapped trays with caps packed separately.

Catalogue No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
FB59543	13	100	13-415	1,000
FB59545	16	100	15-415	1,000
FB59547	16	125	15-415	1,000
FB59549	16	150	15-415	1,000
FB59551	20	125	18-415	500
FB59553	20	150	18-415	500



## Test tubes, borosilicate glass, round bottom, screw thread without caps, disposable 40

**Fisherbrand®**

Ideal for tissue culture, bacteriology, clinical chemistry, blood typing and cross matching procedures.

Available with or without marking spot. Tubes packed in shrink wrapped trays.

With marking spot

Catalogue No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
FB59555	13	100	13-415	1,000
FB59557	16	100	15-415	1,000
FB59559	16	125	15-415	1,000
FB59561	16	150	15-415	1,000
FB59563	20	125	18-415	500
FB59565	20	150	18-415	500

Without marking spot

Catalogue No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
FB59567	13	100	13-415	1,000
FB59571	16	150	15-415	1,000
FB59575	20	150	18-415	500

For caps, see next page.



## Test tubes, borosilicate glass, light walled, rimless 40

**Fisherbrand®**

Premium quality tubes with sturdy, uniform bottoms and consistent lengths. The 10mm x 75mm and 12mm x 75mm sizes are suitable for cell washing procedures. The 6mm x 50mm size is often referred to as a 'Durham tube'. Approximate wall thickness 0.6mm. FB59521 is manufactured from flint glass.

Catalogue No	Diameter [external], mm	Length, mm	Pack qty
FB59521	6	50	1,000
FB59523	10	75	1,000
FB59525	12	75	1,000
FB59527	13	100	1,000
FB59529	15	85	1,000
FB59531	16	100	1,000
FB59533	16	125	1,000
FB59535	16	150	1,000
FB59537	18	150	500
FB59539	20	150	500
FB59541	25	150	500

For caps, see next page.

