


## Tubing, I/P

Pump tubing cross sections



	I/P 73	I/P 82	I/P 70	I/P 88	I/P 89
Inside diameter (nominal), mm	9.5	12.7	9.5	12.7	15.88
Hose barb size, mm	9.5	12.7	9.5	12.7	15.88
Flow range (approximate) <sup>†</sup> with 1 to 650 rpm drive, L/min	0.01 to 0.8	0.02 to 13	0.01 to 8	0.02 to 13	0.03 to 17

<sup>†</sup>Determined at Opsi at inlet, 0.5psi at outlet, water temperature 22°C

CP Tubing ordering information - lengths vary, see our website or contact Customer Service

Tubing type	I/P 26	I/P 73	I/P 82	I/P 70	I/P 88	I/P 89
Silicone (platinum)	PMM-550-510L	PMM-555-510E	PMM-560-510E	PMM-566-510X	N/A	PMM-576-510H
Silicone (peroxide)	PMM-550-515B	PMM-555-515R	PMM-560-515R	PMM-565-515K	PMM-570-515K	PMM-575-515D
BioPharm silicone (platinum)	PMM-550-520Y	PMM-555-520B	PMM-560-520B	PMM-566-520L	PMM-571-520L	PMM-576-520E
BioPharm Plus silicone (platinum)	PMM-550-525V	PMM-555-525X	PMM-560-525X	PMM-566-525B	PMM-571-525B	PMM-576-525R
C-FLEX® (50 A)	PMM-550-530F	PMM-555-530V	PMM-560-530V	PMM-565-530X	PMM-570-530X	PMM-575-530H
PharMed® (65)	PMM-550-500X	PMM-555-500H	PMM-560-500H	PMM-565-500A	PMM-570-500A	PMM-575-500Q
PharmaPure®	PMM-550-502K	PMM-555-502D	PMM-560-502D	N/A	PMM-570-502T	N/A
Sta-Pure®	PMM-550-503Y	PMM-555-503B	PMM-560-503B	N/A	N/A	N/A
Chem-Sure®	PMM-550-504G	PMM-555-504W	PMM-560-504W	N/A	N/A	N/A
BioLine™	N/A	N/A	N/A	N/A	N/A	N/A
Tygon® LFL	PMM-550-505E	PMM-555-505U	PMM-560-505U	PMM-565-505N	PMM-570-505N	PMM-575-505G
Tygon® Food (B-44-4X)	PMM-550-540C	PMM-555-540S	PMM-560-540S	N/A	N/A	N/A
Tygon Lab® (R 3603)	PMM-550-545P	PMM-555-545Y	PMM-560-545Y	PMM-570-545B	N/A	PMM-550-558G
Tygon® Fuel & Lubricant (F-4040-A)	PMM-550-555M	PMM-555-555F	PMM-560-555F	N/A	N/A	N/A
Tygon® Chemical (2001)	PMM-555-561K	PMM-560-562Y	PMM-550-565J	N/A	N/A	N/A
Norprene® (A60 G)	PMM-550-550W	PMM-555-550P	PMM-560-550P	PMM-565-540L	PMM-570-540L	PMM-575-540E
Norprene® food (A60 F)	PMM-550-536Q	PMM-555-536J	PMM-560-536J	PMM-565-535E	PMM-570-535E	PMM-575-535U
Viton®	PMM-550-560T	PMM-555-560M	N/A	N/A	N/A	N/A
FDA Viton®	PMM-550-570Q	PMM-555-570J	N/A	N/A	N/A	N/A