



Vials, sample, clear and amber, Type 1 Class A borosilicate glass, with fitted caps

Fisherbrand®

46

- Supplied in partitioned trays
- Standard vials have fitted black phenolic caps with either polycone liner or Teflon® faced white rubber liner
- Also available: pre-cleaned EPA (Environment Protection Agency) vials for VOA (volatile organic analysis), with or without certification, having fitted white polypropylene open top caps with Teflon® faced silicone septa

Clear sample vials with fitted black phenolic polycone lined caps

Catalogue No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
FB72901	4	1	15	48	13-425	144
FB72902	8	2	17	63	15-425	144
FB72903	12	3	19	68	15-425	144
FB72904	16	4	21	73	18-400	144
FB72907	20	5	28	60	24-400	72
FB72905	24	6	23	88	20-400	144
FB72908	25	6.25	28	73	24-400	72
FB72906	30	8	25	94	22-400	144

Clear sample vials with fitted black phenolic Teflon® faced rubber lined caps

Catalogue No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
FB72911	2	0.5	12	38	8-425	288
FB72912	4	1	15	48	13-425	144
FB72913	8	2	17	63	15-425	144
FB72914	12	3	19	68	15-425	144
FB72915	16	4	21	73	18-400	144
FB72918	20	5	28	60	24-400	72
FB72916	24	6	23	88	20-400	144
FB72919	25	6.25	28	73	24-400	72
FB72917	30	8	25	94	22-400	144

Clear short sample vials with fitted black phenolic Teflon® faced rubber lined caps

Catalogue No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
FB72921	2	0.5	15	31	13-425	200
FB72922	4	1	17	41	15-425	200
FB72923	6	1.5	19	43	15-425	200

Amber sample vials with fitted black phenolic Teflon® faced rubber lined caps

Catalogue No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
FB72925	2	0.5	12	38	8-425	288
FB72926	4	1	15	48	13-425	144
FB72927	8	2	17	63	15-425	144
FB72928	20	5	28	60	24-400	72
FB72929	40	10	28	98	24-400	72

Pre-cleaned EPA (Environmental Protection Agency) vials for VOA in both clear and amber, with fitted white polypropylene open top cap and fitted Teflon® faced silicone septa. All EPA vials have an outside diameter of 28mm, 24-400 GPI thread finish.

Without certification

Catalogue No	Capacity, mL	Capacity, drams	Colour	Height, mm	Pack qty
FB72931	20	5	Clear	60	144
FB72932	20	5	Amber	60	144
FB72933	40	10	Clear	98	144
FB72934	40	10	Amber	98	144

With certification

Catalogue No	Capacity, mL	Capacity, drams	Colour	Height, mm	Pack qty
FB72935	20	5	Clear	60	144
FB72936	40	10	Clear	98	144
FB72937	40	10	Amber	98	144