



Vials, borosilicate glass, large, clear, in Lab File™

03 WHEATON Science International

For storage of reagents and standards. Excellent for small derivatisation reactions. Especially suited where stability of the glass is of prime importance. Packed in convenient Lab File™ cases. Manufactured from Wheaton-33® low extractable borosilicate glass (meets ASTM Type 1 Class A and USP Type 1 standards) to protect against changes in pH and maintain purity of contents.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-835-010N	224831	20	27.8	57	24	400	200
VGA-835-020K	224832	25	27.8	70	24	400	200
VGA-835-030H	224833	40	27.8	95	24	400	200



Vials, borosilicate glass, large, with cap, in Lab File™

03 WHEATON Science International



Specially suited for research studies where stability of the glass is of prime importance. Wheaton-33® low extractable borosilicate glass (meets ASTM Type 1 Class A and USP Type 1 standards). Vials can be autoclaved. Caps have white (styrene-butadiene) rubber liner. Packed in convenient Lab File™ cases.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-840-010N	225288	20	27.8	60	24	400	72
VGA-840-020K	225289	25	27.8	73	24	400	72
VGA-840-030H	225290	40	27.8	98	24	400	72



Vials, borosilicate glass, Shorty Vials® in Lab File™

03 WHEATON Science International

Low profiles for greater stability and less dead volume above sample. Packed in shrink-wrapped trays in convenient Lab File™ cases. Wheaton-33® low extractable borosilicate glass (meets ASTM Type 1 Class A and USP Type 1 standards). Caps not included. See page 1,469.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-845-010G	224821	2	15	28	13	425	200
VGA-845-020D	224822	4	17	38	15	425	200
VGA-845-030A	224823	6	19	40	15	425	200



Vials, storage, with vials and rubber lined caps, Vial File™

03 WHEATON Science International

Vials with white (styrene-butadiene) rubber lined caps and stored in Vial File™.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-850-010G	224891	2	12	38	8	425	60
VGA-850-020D	224892	4	15	48	13	425	40