

Vials, storage, with Teflon® lined caps, Vial File™

03



Open top caps are lined with Teflon® faced silicone rubber. Stored in Vial File™.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-850-030A	224950	1.5	12	35	8	425	60
VGA-850-040U	224952	4	15	48	13	425	40



Vials, V, screw top with Teflon® lined solid top caps

03



Solid top black phenolic screw cap with Teflon® faced 14B white rubber liner.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-855-021R	W986254NG	1.0	13	44	13	425	12
VGA-855-011U	W986253NG	0.3	13	35	13	425	12
VGA-855-030Q	986256	2.0	17	61	15	415	12
VGA-855-041L	W986257NG	3.0	20	50	20	400	12
VGA-855-050K	986259	5.0	20	61	20	400	12



Vials, V, screw top, graduated, with Teflon® lined open top caps

03



Graduated vial with open top black phenolic screw cap and Teflon® faced silicone liner.

Catalogue No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
VGA-860-010W	986293	0.3	13	32	13	425	12
VGA-860-020T	986294	1.0	13	41	13	425	12
VGA-860-030Q	986298	2.0	20	41	20	400	12
VGA-860-040N	986297	3.0	20	47	20	400	12
VGA-860-050K	986299	5.0	20	62	20	400	12



Vials, glass, liquid scintillation

03



The original scintillation vial, 20mL capacity, manufactured from Wheaton '180' low potassium glass under critical manufacturing tolerances. Background counts are consistent and low whilst UV transmission is high. Tight fitting white screw caps supplied. Dimensions 28mm x 61mm (with cap on).

Catalogue No	Alt. No	Closure size, mm	Closure	Cap style	Inner pack qty	Pack qty
VGA-870-010P	986541	22	Polypropylene	Metal foil	100	500
VGA-870-020M	986542	22	Urea	Metal foil	100	500
VGA-870-030J	986546	22	Urea	Poly-seal cone	100	500
VGA-870-040G	986548	22	Urea	Polyethylene disc	100	500
VGA-870-050D	986561	24	Polypropylene	Metal foil	100	500
VGA-870-060A	986568	24	Urea	Polyethylene	100	500



Vials, high density polyethylene, liquid scintillation

03



Wheaton polyethylene vials feature lightweight walls for increased counting efficiency.

They are premium vials, engineered and produced to precise tolerances to avoid jamming. Dimensions 27mm x 61mm (with cap on). Screw cap size 22mm. Capacity 20mL.

Catalogue No	Alt. No	Closure	Cap style	Inner pack qty	Pack qty
VGA-875-010Y	986701	Polypropylene	Metal foil	100	500
VGA-875-020F	986702	Urea	Metal foil	100	500
VGA-875-030C	986704	Polyethylene	Linerless	100	500

