

## Washing machine/disinfector, with HEPA filtered glassware drying unit, Model G 7883 CD

04



Catalogue No	Alt. No
WCK-205-020L	G 7883 CD

### Accessories and spares

Catalogue No	Description
WCK-200-505R	For liquid cleaning agents
WCK-200-510B	Dispenser (DOS G 10) for liquid rinsing agents

### Cleaning agents

Catalogue No	Description	Quantity
FB70173	Cleaning detergent	1L
FB70177	Cleaning detergent	5L
D/0029/21	Decon acid rinse	5L
D/0026/65	Dri-Decon® powder detergent	5L

### Baskets, racks, inserts and other accessories (for G 7804, G 7883 and G 7883 CD)

#### Upper baskets

Catalogue No	Alt. No	Description
WCK-205-500U	O 188	Open front, clearance 165mm, with spray arm
WCK-205-502Q	O 190	Open front, clearance 215mm, with spray arm
WCK-205-504M	O 187	With injector unit, clearance 170mm, for narrow necked glassware (G 7883 only)
WCK-205-505K	O 175	With injector unit, clearance 170mm, for narrow necked glassware (G 7883 CD only)
WCK-205-506Y	O 184	With injector unit, clearance 190mm, for narrow necked glassware (G 7883 only)

#### Lower baskets

Catalogue No	Alt. No	Description
WCK-205-508E	U 874	Open front, clearance 270mm (when used with O 188), 220mm (when used with O 190)
WCK-205-510R	U 175	With injector unit, clearance 170mm, for narrow necked glassware, 33 jets, (for use with O 175 and O 187 upper baskets only)
WCK-205-512N	U 184	With injector unit, clearance 170mm, for tubes and vials, 96 jets, (for use with O 175 and O 184 upper baskets only)

#### Inserts

##### For test tubes

Catalogue No	Alt. No	Description
WCK-206-500X	E 103	For tubes up to 107mm in length, for upper and lower baskets
WCK-206-502K	E 104	For tubes up to 105mm in length, for upper and lower baskets
WCK-206-504G	E 105	For tubes up to 165mm in length, for lower baskets
WCK-206-506C	E 139	For tubes up to 200mm in length, for lower baskets
WCK-206-508V	E 149	For 80 tubes up to 105mm in length, including lid, upper and lower baskets

##### For Erlenmeyer flasks/beakers

Catalogue No	Alt. No	Description
WCK-206-510L	AK 12	Half basket insert, for upper and lower baskets
WCK-206-512H	E 106	With 28 spring clips (10 x 175mm height, 18 x 105mm height), for lower baskets
WCK-206-514D	E 106/1	With 28 spring clips (28 x 105mm height), for upper and lower baskets
WCK-206-516W	E 106/2	With 15 spring clips (15 x 175mm height), for lower baskets

##### For Erlenmeyer flasks/round beakers

Catalogue No	Alt. No	Description
WCK-206-518S	E 109	For 21 beakers to 250mL, for lower baskets
WCK-206-520Y	E 110	For 10 beakers from 250mL to 600mL, for lower baskets
WCK-206-522E	E 111	For 8 beakers from 600mL to 1,000mL, for lower baskets
WCK-206-524A	E 144	for 18 beakers up to 250mL, for upper and lower baskets

entry continued

