

Test sticks, Quantofix®

26



- Rapid: most tests can be carried out in 10 to 120 seconds
- Convenient: ready to use kits, pre-calibrated, contain all necessary equipment and reagents
- Precise: colour charts adjusted and checked using certified standard solutions
- Shelf life: 2.5 years, store below 30°C



| Catalogue No | Alt. No | Description | Graduations | Pack type | Pack qty |
|--------------|---------|-----------------|--|---------------------|----------|
| WAT-700-010F | 91307 | Aluminium | 0, 5, 20, 50, 200, 500 | Sticks and reagents | 100 |
| WAT-700-050Q | 91315 | Ammonium | 0, 10, 25, 50, 100, 200, 400 | Sticks and reagents | 100 |
| WAT-700-700D | 91332 | Arsenic 50 | 0, 0.05, 0.1, 0.5, 1, 1.7, 3 | Sticks and reagents | 100 |
| WAT-700-710A | 91334 | Arsenic 10 | 0, 0.01, 0.025, 0.05, 0.1, 0.5 | Sticks and reagents | 100 |
| WAT-700-060N | 91314 | Ascorbic acid | 0, 50, 100, 200, 300, 500, 1,000, 2,000 | Sticks | 100 |
| WAT-700-070K | 91324 | Calcium | 0, 10, 25, 50, 100 | Sticks and reagents | 60 |
| WAT-700-080H | 91323 | Carbonate | 0, 3, 6, 10, 15, 20d | Sticks | 100 |
| WAT-700-090E | 91321 | Chloride | 0, 500, 1,000, 1,500, 2,000, >3,000 | Sticks | 100 |
| WAT-700-150M | 91301 | Chromate | 0, 3, 10, 30, 100 | Sticks and reagents | 100 |
| WAT-700-160J | 91303 | Cobalt | 0, 10, 25, 50, 100, 250, 500, 1,000 | Sticks | 100 |
| WAT-700-170G | 91304 | Copper | 0, 10, 30, 100, 300 | Sticks | 100 |
| WAT-700-200A | 91318 | Cyanide | 0, 1, 3, 10, 30 | Sticks and reagents | 100 |
| WAT-700-210U | 91335 | EDTA | 0, 100, 200, 300, 400 | Sticks | 100 |
| WAT-700-220R | 91328 | Formaldehyde | 0, 10, 20, 40, 60, 100, 200 | Sticks and reagents | 100 |
| WAT-700-305J | 91308 | Iron 100 | 0, 2, 5, 10, 25, 50, 100 | Sticks and reagents | 100 |
| WAT-700-300T | 91302 | Iron 1000 | 0, 5, 20, 50, 100, 250, 500, 1,000 | Sticks and reagents | 100 |
| WAT-700-320N | 91336 | LubriCheck | 0, 15, 50, 75, 130, 200mmol/L IKOH | Sticks | 100 |
| WAT-700-330K | 91325 | Molybdenum | 0, 5, 20, 50, 100, 250 | Sticks and reagents | 100 |
| WAT-700-350E | 91305 | Nickel | 0, 10, 25, 50, 100, 250, 500, 1,000 | Sticks | 100 |
| WAT-700-400P | 91313 | Nitrite/nitrate | 0, 10, 25, 50, 100, 250, 500 for NO ₃ 0, 1, 5, 10, 20, 40, 80 for NO ₂ | Sticks | 100 |
| WAT-700-450A | 91311 | Nitrite | 0, 1, 5, 10, 20, 40, 80 | Sticks | 100 |
| WAT-700-460U | 91338 | Nitrite pH | 0, 1, 5, 10, 20, 40, 80 for NO ₂ pH 6.0, 6.4, 6.7, 7.0, 7.3, 7.6, 7.9, 8.2, 8.4, 8.6, 8.8, 9.0, 9.3, 9.6 | Sticks | 100 |
| WAT-700-470R | 91322 | Nitrite 3000 | 0, 0.1, 0.3, 0.6, 1, 2, 3g/L | Sticks | 100 |
| WAT-700-480X | 91333 | Peroxide 1000 | 0, 50, 150, 300, 500, 800, 1,000 | Sticks | 100 |
| WAT-700-485E | 91312 | Peroxide 100 | 0, 1, 3, 10, 30, 100 | Sticks | 100 |
| WAT-700-490L | 91319 | Peroxide 25 | 0, 0.5, 2, 5, 10, 25 | Sticks | 100 |
| WAT-700-495B | 91320 | Phosphate | 0, 3, 10, 25, 50, 100 PO ₄ | Sticks and reagents | 100 |
| WAT-700-500L | 91316 | Potassium | 0, 200, 400, 700, 1,000, 1,500 | Sticks and reagents | 100 |
| WAT-700-510Y | 91337 | Quat | 0, 10, 25, 50, 100, 250, 500, 1,000 | Sticks | 100 |
| WAT-700-610E | 91329 | Sulfate | <200, >400, >800, >1,200, >1,600 | Sticks | 100 |
| WAT-700-600H | 91306 | Sulfite | 0, 10, 25, 50, 100, 250, 500, 1,000 | Sticks | 100 |
| WAT-700-550T | 91309 | Tin | 0, 10, 25, 50, 100, 250, 500 | Sticks | 100 |
| WAT-700-620B | 91310 | Zinc | 0, 2, 5, 10, 25, 30, 100 | Sticks and reagents | 100 |