



Thermo Scientific
Nunc Serological Pipettes

research with confidence



Thermo Scientific Nunc Serological Pipettes

Our newly expanded portfolio of serological pipettes features a full range of size and packaging options – offering the quality, accuracy and performance today's scientists require during every stage of cell culture research.

Cell isolation
and preparation



Characterization
and Analysis



We are dedicated to meeting the rigorous demands of your research, at every stage of cell culture. Our new collection of serological pipettes offers the quality assurance, size variability and packaging options you need to advance your research with complete confidence.

Growth
and passage



Accurate, disposable plastic pipettes plugged and sterilized

- Sterility assurance level (SAL) of 10^{-6}
- 100% virgin polystyrene for maximum clarity
- Certified non-pyrogenic
- Bold clear graduations with easy-to-read scale markings
- Non-dripping tip for pipette precision
- Color-coded packaging for ease in sorting and selecting the correct size

Packaging options that accommodate every lab!

- Individually wrapped; paper/plastic or plastic film; peel or pop through for easy opening
- Bulk packed for frequent use; reduces packaging waste
- Multiple sizes for use in a wide range of applications

Storage

Harvest

Serological Portfolio Highlights

- Choose multiple sizes, ranging from 1 mL to 50 mL, all with non-dripping tips and clearly marked graduations for precision
- Pick bulk packed options for higher volume use with less packaging waste
- Select individually wrapped pipettes in either paper/plastic peel or plastic film
- See the difference with new ultra clear, 100% virgin polystyrene



Thermo Scientific Nunc Serological Pipettes

Ordering Information

Cat No	Manu. No	Description	Unit per bag	Unit per case
Individually wrapped, paper peel packaging, plugged				
13450287	170353	1 ml pipette	100	1000
13470287	170354	2 ml pipette	100	500
13430287	170355	5 ml pipette	50	200
13440287	170356	10 ml pipette	50	200
13460287	170357	25 ml pipette	50	200
13420297	170358	50 ml pipette	25	100
Shortie, Individually wrapped paper/plastic peel packaging, plugged				
15101539	170360	5 ml shortie pipette	50	200
15191529	170361	10 ml shortie pipette	50	200
Individually wrapped, plastic film packaging, plugged				
13480287	170364	1 ml pipette	100	1000
13490287	170365	2 ml pipette	100	500
13400297	170366	5 ml pipette	50	200
13410297	170367	10 ml pipette	50	200
13430297	170368	25 ml pipette	50	200
13440297	170369	50 ml pipette	25	100
Pipettes, Bulk Packed, Plugged				
13450297	170371	1 ml pipette	50	1000
13460297	170372	2 ml pipette	50	500
13470297	170373	5 ml pipette	25	500
13480297	170374	10 ml pipette	25	500
13490297	170375	25 ml pipette	25	200
13400307	170376	50 ml pipette	25	100



Thermo Scientific S1 Pipet Filler

A lightweight, cordless pipet filler that allows longer, fatigue-free pipetting—an excellent companion to our serological pipettes.

For more information, go to
www.thermoscientific.com/s1pipetfiller

© 2014 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 **Belgium:** +32 (0)56 260 260 **Denmark:** +45 70 27 99 20
Germany: +49 (0)180 5258221 **Ireland:** +353 (0)1 885 5854 **Italy:** +39 02 950 59 478
Finland: +358 (0)9 8027 6280 **France:** +33 (0)3 88 67 14 14 **Netherlands:** +31 (0)20 487 70 00
Norway: +47 22 95 59 59 **Portugal:** +351 21 425 33 50 **Spain:** +34 902 239 303
Sweden: +41 31 352 32 00 **Switzerland:** +41 (0)56 618 41 11 **UK:** +44 (0)1509 555 500

 **Fisher
Scientific**
Part of Thermo Fisher Scientific