

Bottles, borosilicate glass, reagent and media, screw neck



PYREX®

Graduated with black phenolic screw cap with PTFE faced sealing disc for chemical inertness.

Resistant to repeated dry heat and steam sterilisation cycles (121°C for 2hr). Autoclavable, do not autoclave with tightly screwed cap.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	O.D., mm	ISO thread GL	Pack qty
11940326	-	1515/01D	25	65	35	20	10
11961598	-	1515/02D	50	80	43	25	10
10105044	-	1515/04D	100	105	53	25	10
10411235	-	1515/06D	250	135	66	30	10
10786885	-	1515/08D	500	175	81	30	10
11337854	-	1515/10D	1,000	215	101	42	10
13115066	-	1515/12D	2,000	245	136	60	10
11347854	-	1515/16D	5,000	324	181	60	1
11357854	-	1515/18D	10,000	395	227	60	1
11367854	-	1515/20D	15,000	460	258	60	1
11377854	-	1515/22D	20,000	502	288	60	1



Black phenolic plastic caps, with PTFE faced sealing disc of ethylene propylene, as described above

Catalogue No	Legacy No	Alt. No	Cap size, mm	For bottle capacity, mL	Pack qty
10399974	-	4503/20	20	25	20
10034794	-	4503/25	25	50, 100	20
12579835	-	4503/30	30	250, 500	20
10181930	-	4503/42	42	1,000	10
10417050	-	4503/60	60	2,000 to 20,000	5



PTFE faced sealing disc for above

Catalogue No	Legacy No	Alt. No	Cap size, mm	Pack qty
10480465	-	4513/20	20	20
10034314	-	4513/25	25	20
10325040	-	4513/30	30	20
10213340	-	4513/42	42	10
11393484	-	4513/60	60	5

Bottles, borosilicate glass, reagent and media, screw neck, ISO 4796, with retrace code



Graduated with blue polypropylene screw cap, ISO thread.

Supplied with interchangeable cap and drip-resistant plastic pouring ring. These bottles are stable, chemically resistant, and can be autoclaved up to 140°C. Do not autoclave bottles with a tightly screwed cap.

For spare caps and pouring rings, see page 175.

Catalogue No	Legacy No	Alt. No	Capacity, mL	ISO thread GL	Height, mm	O.D., mm	Pack qty
13117384	-	218011453	25	25	74	36	10
12307178	-	218011753	50	32	91	46	10
12317178	-	218012458	100	45	105	56	10
13127216	-	218012955	150	45	115	62	10
12327178	-	218013651	250	45	143	70	10
12337178	-	218014459	500	45	181	86	10
13137216	-	218015155	750	45	208	95	10
12347178	-	218015455	1,000	45	230	101	10
12357178	-	1018630	2,000	45	267	136	1
13147216	-	218016957	3,500	45	300	160	1
13197484	-	218017353	5,000	45	340	186	1
13134704	-	218018658	10,000	45	415	234	1
13197843	-	218019159	20,000	45	510	300	1

