

Cuvettes, disposable, UV, micro, UV-Cuvette



- Disposable plastic cuvette, no carry over of samples, reduced risk of contamination
- UV and Vis transparent between 220nm and 900nm
- Minimum sample 70µL
- Ideal for determination of proteins, DNA, RNA in the UV range
- Fits into most commercial photometers
- Two centre heights available: 8.5mm and 15mm
- Suitable for freezing down to -20°C when using caps
- Dimensions of window (w x h) 2mm x 3.5mm
- Max 70°C
- Individually wrapped cuvettes are free of DNase, DNA and RNase

UV Cuvette micro cuvettes are made from UV transparent material ideally suited for applications in biology and biochemistry. Replace expensive and fragile glass or quartz glass cuvettes. A round cap is available for tight sealing. 10730445 and 10386712 are individually wrapped.



Catalogue No	Legacy No	Alt. No	Capacity, µL	Height, mm	Pack qty
10290272	-	759200	70 to 850	8.5 (centre)	100
10188561	-	759210	70 to 850	8.5 (centre)	500
10730445	-	759215	70 to 850, individually wrapped	8.5 (centre)	100
10200282	-	759220	70 to 550	15 (centre)	100
10198561	-	759230	70 to 550	15 (centre)	500
10386712	-	759235	70 to 550, individually wrapped	15 (centre)	100

Accessory caps, PE round cover

Catalogue No	Legacy No	Alt. No	Colour	Pack qty
10086751	-	759240	Blue	100
10159441	-	759241	Yellow	100
10476221	-	759242	Green	100
10139441	-	759243	Orange	100

Cuvettes, disposable, semi-micro and macro



- Path length 10mm
- PS and PMMA grouped by mould cavity number to ensure lowest variation of extinction coefficient

UV cuvettes have a very high chemical resistance. They are ideal for the determination of proteins, DNA and RNA. They fit virtually all commercial photometers for standard cuvettes.

Technical Specification - Specific

	UV	PS	PMMA
Wavelength range, nm	220 to 900	340 to 900	300 to 900

Catalogue No	Legacy No	Alt. No	Type	Description	Capacity, mL	Pack qty
10566581	-	759170	UV	Semi-micro	2.5 to 4.5	100
10046731	-	759150	UV	Macro	1.5 to 3.0	100
12741735	-	759071D	PS	Macro	2.5 to 4.5	100
12328875	-	759076D	PS	Semi-micro	1.5 to 3.0	100
12761735	-	759081D	PMMA	Macro	2.5 to 4.5	100

