

# Plastic Transfer Pipettes

- LDPE, Low density polyethylene
- Transparent
- Graduated or non-graduated
- Exists in sterile version
- Various packaging




Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu\text{L}$	Drop per ml	Quantity
13469118	Transfer Pipette PE, 1mL	No	104	33	30	400



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu\text{L}$	Drop per ml	Quantity
13499108	Transfer Pipette PE, 1mL, Graduated	No	150	33	30	500



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu\text{L}$	Drop per ml	Quantity
13439118	Transfer Pipette PE, 1mL, Graduated	No	150	33	30	500
13489108	Transfer Pipette PE, 1mL, Sterile, Graduated	Yes	150	33	30	500
13419118	Transfer Pipette PE, 1mL, Sterile, Graduated, inner pack of 10	Yes	150	33	30	500
13429118	Transfer Pipette PE, 1mL, Sterile, Graduated, inner pack of 20	Yes	150	33	30	500



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu\text{L}$	Drop per ml	Quantity
13439108	Transfer Pipette PE, 3mL, Graduated	No	155	40	25	500
13469108	Transfer Pipette PE, 3mL, Sterile, Graduated	Yes	155	40	25	500
13479108	Transfer Pipette PE, 3mL, Sterile, Graduated, inner pack of 10	Yes	155	40	25	500
13409118	Transfer Pipette PE, 3mL, Sterile, Graduated, inner pack of 20	Yes	155	40	25	500

**Terms and Conditions** Please note that all sales are subject to our standard terms and conditions. Full details can be found on our website [www.eu.fishersci.com](http://www.eu.fishersci.com)



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu$ L	Drop per ml	Quantity
13459118	Transfer Pipette PE, 4mL, thin stem	No	150	33	30	500



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu$ L	Drop per ml	Quantity
13459108	Transfer Pipette PE, 4mL	No	150	33	50	500



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu$ L	Drop per ml	Quantity
13449108	Transfer Pipette PE, 7mL, extra long	No	300	50	20	100



Cat. No	Description	Sterile	Length, mm	Drop volume, $\mu$ L	Drop per ml	Quantity
13449118	Transfer Pipette PE, 10mL, jumbo	No	170	56	18	200

**Terms and Conditions** Please note that all sales are subject to our standard terms and conditions. Full details can be found on our website [www.eu.fishersci.com](http://www.eu.fishersci.com)