



Gloves, latex, Class 1000, ISO Class 6,
hand specific PPE Category III,
AccuTech® Gammex®

Ansell

S3



AccuTech® Gammex® sterile gloves are made from specific formulation of the latex that provides increased flexibility and elasticity.

The optimised thickness of the latex, combined with the bisque finish, provides outstanding fingertip sensitivity plus increased handling precision and dexterity. This makes these gloves ideal for intricate handling operations. Length 300mm. Pack of 200 pairs.

Catalogue No	Size	Pack type	Pack qty
SAR-467-010N	6	Pair	200
SAR-467-030H	6.5	Pair	200
SAR-467-050B	7	Pair	200
SAR-467-070S	7.5	Pair	200
SAR-467-090M	8	Pair	200
SAR-467-110J	8.5	Pair	200
SAR-467-130D	9	Pair	200

Gloves, latex, powder free, surgeons, Profile 3000™

04



Hand specific glove manufactured to exceed requirements of European Medical Glove Standard EN 455 (1 & 2). Conforms to ISO 9002 and EN 46002. Gamma irradiated, bar-coded and lot numbered for traceability. Pairs individually packaged in a sealed pouch, for use in controlled environments. Anatomically shaped and textured for excellent grip, particularly when wet. Length 305mm.

Catalogue No	Size	Pack type	Pack qty
CMK-140-010Q	9.5	Pair	50
CMK-140-020N	6	Pair	50
CMK-140-030K	6.5	Pair	50
CMK-140-040H	7	Pair	50
CMK-140-050E	7.5	Pair	50
CMK-140-060B	8	Pair	50
CMK-140-070V	8.5	Pair	50
CMK-140-080S	9	Pair	50
CMK-140-090P	10	Pair	50

Gloves, latex, powder free, BioClean Class 100 ISO
Class 5 cleanroom compatible

04



Hand specific glove made to the same specification as CMK-100-range. Particularly suitable for use in environments where fluid contamination may create a serious risk. Length 610mm.

Catalogue No	Size	Pack type	Pack qty
CMK-115-020E	6	Pair	25
CMK-115-040V	7	Pair	25
CMK-115-050S	7.5	Pair	25
CMK-115-060P	8	Pair	25
CMK-115-070M	8.5	Pair	25
CMK-115-080J	9	Pair	25

Gloves, latex, powder free, BioClean Class 100
ISO Class 5 cleanroom compatible ambidextrous

04



Packed in sealed plastic inner bags and further wrapped in an outer sealed bag for cleanroom compatibility. Textured for improved grip. Particulates and extractables are minimised by extended washing. Economical and extremely low in latex allergens. Manufactured to exceed the requirements of the European Medical Glove Standard EN 455 (Parts 1 & 2). Ambidextrous. Length 305mm.

Catalogue No	Size	Pack qty
CMK-165-030M	Small	100
CMK-165-040J	Medium	100
CMK-165-050G	Large	100
CMK-165-060D	XL	100