



Gloves, co-polymer, Profile P-Zero™, cleanroom,  
ISO Class 5, Class 1 Medical

S4

**S**



- EN 455 (parts 1 & 2)
- AQL 1.5
- Compatible for use in ISO Class 5 cleanrooms
- Double bagged
- Sterile
- 300mm (12 inch) length for extra protection
- Hand-specific

Catalogue No	Size	Pack type	Pack qty
CMK-450-020F	5.5	Pair	200
CMK-450-030C	6.0	Pair	200
CMK-450-040W	6.5	Pair	200
CMK-450-050T	7.0	Pair	200
CMK-450-060Q	7.5	Pair	200
CMK-450-070N	8.0	Pair	200
CMK-450-080K	8.5	Pair	200
CMK-450-090H	9.0	Pair	200

Gloves, nitrile, cleanroom,  
HypoClean Critical™

04

**Kimberly-Clark**  
Professional

**S**

Designed to exceed the stringent requirements of state of the art critical environment facilities, this new technology combines industry leading contamination control performance with comfort and durability to create a new benchmark in cleanroom gloves. Ultra-clean performance is assured each time you don the HypoClean Critical™ glove. The proprietary compound, customised processes and rigid statistical process control ensure optimum consistency. Vertically integrated production and testing in ISO 9002 certified facilities enable production of a glove with the lowest levels of particle and extractable contamination, eliminating the for incoming inspection. Length 305mm.



Catalogue No	Size	Pack qty
CMK-350-020N	XS	20
CMK-350-040H	Small	20
CMK-350-060B	Medium	20
CMK-350-080S	Large	20
CMK-350-100P	XL	20

Gloves, nitrile, N-Plus™, cleanroom ISO Class 5,  
Class 1 Medical

S4

**S**



- EN 455 (parts 1 & 2)
- AQL 1.5
- Compatible for use in ISO Class 5 cleanrooms
- Double bagged
- Sterile
- 400mm (16 inch) length for extra protection
- Hand-specific

Catalogue No	Size	Pack type	Pack qty
CMK-210-010G	6.0	Pair	200
CMK-210-020D	6.5	Pair	200
CMK-210-030A	7.0	Pair	200
CMK-210-040U	7.5	Pair	200
CMK-210-050R	8.0	Pair	200
CMK-210-060X	8.5	Pair	200
CMK-210-070L	9.0	Pair	200
CMK-210-080Y	10.0	Pair	200

Gloves, nitrile, non-sterile,  
cleanroom, HypoClean Critical™

04

**Kimberly-Clark**  
Professional

Superior cleanliness, low extractable levels, ideal electrostatic dissipative properties, and superb abrasion resistance make this glove the ideal choice for controlled environments. Prolonged wear, comfort, tactile sensitivity and easy donning are just some of the unique characteristics of the virtually odourless nitrile compound. Beaded cuff and textured fingertips. Ambidextrous. Length 305mm.



Catalogue No	Size	Pack qty
CMK-300-020C	XS	100
CMK-300-040T	Small	100
CMK-300-060N	Medium	100
CMK-300-080H	Large	100
CMK-300-100E	XL	100