

Precision Portable balances

Scout® STX Series



Setting New Standards in Laboratory & Industrial Weighing – The Next Generation of Scout Portable Balances

- Navigate Quickly & Efficiently Using the Informative Color Touchscreen Display
- Fast Weighing Speed and High Resolution Deliver Repeatable and Reliable Results
- Superior Overload Protection & Stackable Storage Drive Overall Durability
- Multiple Weighing Modes and Languages for Simplified Weighing

Application

Weighing, Percent Weighing, Parts Counting, Check Weighing, Dynamic/Animal Weighing, Totalization, Density Determination, Display Hold, Mole Weighing

Display

Full-color Touchscreen WQVGA Graphic LCD (4.3 in. diagonal), QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data

Operation

AC Adapter (Included), batteries (not included)

Communication

Optional RS232, USB host, USB device, Bluetooth or Ethernet

Construction

ABS housing, stainless steel pan, replaceable in-use cover

Design Features

Fast stabilization time, up-front level indicator, up to 9 operating languages, save-to-USB function, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Model	Capacity	Readability (d)	Platform	Test weight	Item №
STX123	120g	0,001g	Ø 93 mm	100g, M1	30253058
STX223	220g	0,001g	Ø 93 mm	200g, M1	30253059
STX222	220g	0,01g	Ø 120 mm	200g, M1	30253060
STX422	420g	0,01g	Ø 120 mm	200g, M1	30253061
STX622	620g	0,01g	Ø 120 mm	300g, M1	30253062
STX1202	1200g	0,01g	170mm x 140mm	1000g, M1	30253063
STX2202	2200g	0,01g	170mm x 140mm	2000g, M1	30253064
STX421	420g	0,1g	Ø 120 mm	200g, M1	30253065
STX621	620g	0,1g	170mm x 140mm	300g, M1	30253066
STX2201	2200g	0,1g	170mm x 140mm	2000g, M1	30253067
STX6201	6200g	0,1g	170mm x 140mm	5000g, M1	30253068
STX8200	8200g	1g	170mm x 140mm	8000g, M1	30253069

Accessories	Item №
RS232 Interface Kit, Scout	30268982
USB Host Interface Kit, Scout	30268983
USB Device Interface Kit, Scout	30268984
Bluetooth Interface Kit, Scout	30268985
Ethernet Interface Kit, Scout	30268986
Stacking Kit, set (6), Scout	30268987
Stacking Kit, set (1), Scout	30268988
Full-featured Impact Printer SF40A	30064202
SF40A Paper roll (57,5mm 2pcs)	12120799
SF40A Ink ribbon cassette	12120798
Auxiliary Display Kit, Scout	30269019
Security device (Laptop lock)	80850043
Top loading kit for Density Determination	30269020
Carrying case, Scout	30269021
In-use cover, Scout	30269022



SF40A IMPACT PRINTER