

SPECIFICATION – CG refrigerated immersion coolers

			C1G	C2G
Temperature range	ambient 20°C	°C	0 to 35°C	-20 to 35°C ¹
Cooling power	@ 20°C	W	300	350 ²
	@ 0°C	W	100 ²	250 ²
	@ -10°C	W	–	150 ²
Overall consumption		W	200	400
Weight		kg	20	24
Dimensions	d/w/h	mm	410/285/225	475/320/255
Flexible pipe	l	mm		925
Coil	Ø/l	mm	77/55	75/105
Electrical supply		V		220-240

¹C1G and C2G immersion coolers will achieve the specified lowest temperatures provided they are used in conjunction with a sufficiently well-insulated bath.

²The nett wattage achievable when used with a circulating bath varies according to tank size and ambient temperature.