



TLC plates, polyester, silica gel 60

- Pre-coated plates 50mm x 200mm, layer thickness 0.25mm
- Contains indicator which fluoresces at 254nm
- Supplied in packs of 50

Catalogue No
TLC-410-610X

Chromatography **Fisherbrand** paper

- Fisherbrand CHR 100 grade of paper may be used in many applications where Whatman No. 1 CHR or Schleicher & Schuell No. 2043a are used
- Fisherbrand CHR 200 grade of paper may be used in many applications where Whatman No. 3MM CHR or Schleicher & Schuell No. 4369 are used
- Supplied in packs of 100

Catalogue No	Dimensions, mm
Fisherbrand CHR 100 FB59501	200 x 200
Fisherbrand CHR 200 FB59515	460 x 570



Aerosol spray bottles, HDPE

- For spraying various liquids, using air as the propellant
- Pressurised by pump action through a unique two way valve
- Can be refilled and repressurised as required, or depressurised for storage or transportation
- Height without cap 165mm, with cap 180mm, diameter 45mm
- Comprises an HDPE bottle and white cap/pump
- Supplied in packs of 5

Catalogue No Capacity, mL
AET-105-Q 200



**Hot air blower, Drying,
General Labware, page 422**

Chromatography paper

- Grade No 1 general purpose paper with a medium flow rate
- For blotting techniques chromatography and electrophoresis
- Available in reels or sheets

Catalogue No	Width, mm
Reels, 100m	
CJF-230-010G	10
CJF-230-020D	20
CJF-230-030A	30
CJF-230-040U	40
CJF-230-050R	50
CJF-230-090F	100
CJF-230-120W	150

Slotted sheets, 12 strips, 15mm wide, pack of 100

Catalogue No	Dimensions, mm
CJF-250-030J	200 x 20
CJF-250-040G	250 x 250
CJF-250-044V	300 x 100
CJF-310-P	213 x 110

