



Funnels, borosilicate glass,
conical, plain


Fisherbrand®

Catalogue No	Top dia., mm	Pack qty
FB56640	55	10
FB56641	75	10
FB56642	100	10
FB56643	150	each
FB56644	200	each



Funnel, Pyrex®
glass, conical,
powder

Catalogue No	Top dia., mm	Stem length, mm	Stem O.D., mm
FPH-700-030U	80	25	18



Funnels, Pyrex® glass, conical,
plain, short stem

Catalogue No	Top dia., mm	Stem length, mm	Stem O.D., mm
FPH-360-005J	35	35	8
FPH-360-010Q	45	45	8
FPH-360-020N	65	55	8
FPH-360-030K	75	75	8
FPH-360-050E	100	100	10



Funnels, polypropylene, conical, plain

Catalogue No	Top dia., mm	Pack qty
FPH-450-010P	40	10
FPH-450-030J	63	10
FPH-450-050D	90	10
FPH-450-070U	115	10
FPH-450-090X	140	5
FPH-450-111J	165	5
FPH-450-130F	200	each



Funnels, PMP, conical

- Internal helicoid ribs to increase the filtration speed of filter paper
- Vertical outer ribs prevent airlock
- Unaffected by corrosive liquids
- Heat resistant to 177°C (350°F)

Catalogue No	Top dia., mm	Stem I.D., mm	Stem O.D., mm	Stem length, mm	Pack qty
FPH-430-030A	51	3	7.5	150	6
FPH-430-050R	70	3	7.5	150	6
FPH-430-070L	70	7	10.0	70	6
FPH-430-090F	100	7	10.0	108	4
FPH-430-110C	140	10	16.0	132	3
FPH-430-130T	196	20	26.0	160	2



Funnel, HDPE,
conical, plain

- With hanging eye for storage
- Ribbed on outside to eliminate back pressure
- Supplied in packs of 5

Catalogue No	Top dia., mm	Vol., mL	Upstand, mm
FPH-470-017K	120	450	15
FPH-485-040B	180	1,500	20
FPH-485-050V	230	2,750	25