



Sharps containers, UN approved, BS7320, UN 3291

NEW

- Sharps containers conforming to BS and UN Standards relating to the transportation of dangerous goods
- Can withstand drop tests at sub zero temperatures
- Visible fill line
- NOTE: Consignor of clinical waste must ensure appropriate packaging is used

Catalogue No	Cap., L	Dims., ext. (w x d x h), mm
SAT-590-011E	0.6	85 x 46 x 195
SAT-590-031V	6	220 x 220 x 230
SAT-590-051P	12	260 x 260 x 325
SAT-590-071J	22	315 x 315 x 385



Disposal, glass

Fisherbrand

- Constructed from corrugated cardboard with integral lid and thick polyethylene bag liner designed for separating glass waste from other types of waste
- Easy to assemble flat pack
- When full the lid is removed and the box is sealed via a rear flap, ready for disposal
- Supplied in packs of 6

Catalogue No	Description	Dims., ext. (w x d x h), mm
FB55141	Floor type	300 x 300 x 680
FB55142	Bench type	200 x 200 x 260



Disposal, glass

- Highly visible receptacle that helps reduce injuries by providing a safe laboratory work area
- Made of sturdy corrugated cardboard, supplied with a polyethylene bag to contain broken glass fragments and spilled liquids.
- When full the internal safety cap is pulled into place and the entire unit is ready for disposal
- Easy to assemble flat pack
- Supplied in packs of 6

Catalogue No	Description	Dimensions, ext. (w x d x h), mm
SAT-625-010R	Floor type	300 x 300 x 680
SAT-625-020X	Bench type	200 x 200 x 260



Hazardous Spill Response kits

- Kits are to be deployed wherever hazardous liquids are in use, for the protection of both personnel and equipment
- In the event of a spill the kit can be used to contain and clean up the hazardous liquid without delay

Contents:

	SAT-670-P	SAT-675-Y
Chemical sorbent sheets (280mm x 330mm)	10	10
Chemical sorbent pillow (230mm x 380mm)	1	5
Chemical sorbent mini booms (1,200mm x 75mm)	-	5
Temporary disposal bag with tie	1	2
Resealable polyethylene bin with lid, 55L	-	1
Wall mounting bracket	-	1
Instruction guide on disposal	1	1

Catalogue No	Cap., L
SAT-670-P	5
SAT-675-Y	26