

Thermometers


Fisherbrand®



Thermometers, white backed, red spirit filled

- General purpose/stirring thermometers with amber graduations
- Reinforced bulb for stirring
- 76mm immersion
- Non-roll fitting

Catalogue No	Range, °C	Subdivisions, °C	Length, mm
FB58547	-10 to +50	0.5	305
FB58549	-10 to +110	1	305
FB58551	-10 to +150	1	305



Thermometers, white backed, blue spirit filled

- General purpose/stirring thermometers with amber graduations
- Available as total immersion or 76mm immersion
- Low toxic environmentally friendly liquid
- Manufactured to traceable standards
- Non-roll top and suspension ring

Catalogue No	Range, °C	Subdivisions, °C	Length, mm
Total immersion			
THL-460-010D	-10 to +210	1	405
THL-460-060L	-20 to +150	1	305
THL-460-120T	-10 to +50	1	155
THL-460-040R	-10 to +50	0.5	305
THL-460-050X	-20 to +110	1	305
THL-460-110W	-10 to +110	1	155
76mm immersion			
THL-460-080F	-10 to +50	0.5	305
THL-460-020A	-10 to +260	1	405
THL-460-090C	-20 to +110	1	305

Thermometers, yellow backed, mercury filled

- Melting point thermometers with black graduations
- Gas filled with expansion chamber
- 30mm immersion
- Non-roll fitting

Catalogue No	Range, °C	Subdivisions, °C	Length, mm
THR-690-030P	0 to +150	0.5	250
THR-690-070D	+140 to 300	0.5	250
THR-690-090U	+290 to 400	0.5	250
THR-690-110R	0 to 360	1	300

Fax: 01509 231893

Telephone: 01509 233344