

NOTE: All timers supplied with an ISO 17025, A2LA Traceable® certificate to NIST (National Institute of Standards and Technology)



**Timer, Big-Foot, Traceable®**

- Big-Foot digits may be read from 9m
- Simultaneous programming of three timing channels
- Functions as a stopwatch with resolution of 1s
- Accuracy 0.01%
- Second super-jumbo display shows time of day

Catalogue No  
FB70220

Spare Part  
BMD-120-022N Lithium manganese  
dioxide batteries 3V  
(CR2025), pack of 5



**Timer, lab, Traceable®**

- Counts down from 9hr, 59min, 59s
- 3 channels, each with distinct alarm
- Accuracy 0.01%
- Dimensions, mm: 75 x 75 x 34
- Supplied with batteries

Catalogue No  
FB70222

Spare Part  
BMT-680-050T Alkaline batteries  
1.5V (LR6), pack of 4



**Timer, 100 hour, Traceable®**

- Times up/down to 99hr, 59min, 59s
- Channels are 4 countdown, 1 clock and 1 stopwatch
- Alarm sounds for 1min, display flashes 'TIME IS UP'
- Accuracy 0.01%
- Dimensions, mm: 70 x 60 x 16
- Supplied with battery

Catalogue No  
FB70224

Spare Part  
BMD-140-011E Silver oxide batteries  
1.55V (SR44),  
pack of 4



**Timer, 3 channel alarm, Traceable®**

- Adjustable alarm volume for quiet or noisy areas
- Program times from 1s to 100hr with memory
- 'Continuous alarm' setting ensures alarm is never missed
- Accuracy 0.001%

Catalogue No  
FB70228

Spare Part  
BMT-680-040W Alkaline batteries  
1.5V (LR03),  
pack of 4



**Timer, 4 channel alarm, Traceable®**

- Easy viewing digits are 19mm high
- Count down/up timing on four separate channels
- Program times from 1s to 100hr with memory
- Displays time-of-day
- Accuracy 0.01%
- Dimensions, mm: 69 x 38 x 12.5

Catalogue No  
FB70230  
Spare Part

BMD-140-011E Silver oxide batteries  
1.55V (SR44),  
pack of 4



**Timer, extra loud, Traceable®**

- May be heard in the noisiest lab or plant environments
- 100dB alarm
- Two timing channels may be used simultaneously
- Attach with magnetic back or use the flip open stand
- Accuracy 0.01%
- Dimensions, mm: 100 x 63 x 12.5

Catalogue No  
FB70232

Spare Part  
BMT-680-040W Alkaline batteries  
1.5V (LR03),  
pack of 4