



**Fisherbrand**<sup>®</sup>

SureOne™ Pipette Tips






*Comprehensive line of non-filter and aerosol barrier pipette tips*



**Special Prices 2009**

# Packaging Options

## Packaging Options

					
	<b>Bulk</b>	<b>Racked</b>	<b>Racked Sterile</b>	<b>Reload</b>	<b>Empty Rack</b>
Pack Format	1000 tip x bag	96 tips x tray 10 trays x pk	96 tips x tray 10 trays x pk	96 tips x tray 10 trays x pk	Inserts not included 10 empty trays x pk
Case Format	10 bags x cs	5 pks x cs	5 pks x cs	5 pks x cs	5 pks x cs









SureOne Standard Sterility Claims - Nonsterile products are certified to be free of RNase, DNase and DNA. E-beam sterile products are certified to be free of RNase, DNase, Pyrogen, Bioburden, PCR Inhibitors, Endotoxin and DNA.

## SureOne Filter Tips

Tip Volume	Packaging	Cat. No.	Alt. No.	Pack Price	
0.1 - 10 µL	Filtered Sterile	FB78098	02-707-439	£42.00	 Micropoint tip, graduated at 2 µL
0.1 - 10 µL	Filtered Sterile	FB78100	02-707-442	£42.00	 Micropoint tip graduated at 2.5 µL extended length
0.1 - 10 µL	Filtered Sterile	FB78102	02-707-456	£42.00	 Micropoint tip, Eppendorf style
2 - 20 µL	Filtered Sterile	FB78104	02-707-435	£42.00	 Micropoint tip, Eppendorf style
2 - 20 µL	Filtered Sterile	FB78106	02-707-432	£42.00	 Beveled Tip
10 - 100 µL	Filtered Sterile	FB78108	02-707-431	£42.00	 Beveled Tip
20 - 200 µL	Filtered Sterile	FB78110	02-707-430	£42.00	 Micropoint tip graduated at 10, 50, 100 µL
100 - 1000 µL	Filtered Sterile	FB78112	02-707-404	£46.00	 Micropoint tip graduated at 100, 200, 500, 1000 µL


# Packaging Options

## SureOne Standard Tips


Tip Volume	Packaging	Cat. No.	Alt. No.	Pack Price	Color	
0.1 - 10 µL	Bulk	FB78000	02-707-454	£8.00	Clear	 Micropoint tip, graduated at 2 µL
	Racked	FB78002	02-707-436	£15.50	Clear	
	Racked Sterile	FB78004	02-707-437	£18.00	Clear	
	Reload	FB78006	02-707-438	£14.85	Clear	
0.1 - 10 µL	Bulk	FB78008	02-707-455	£9.00	Clear	 Micropoint tip graduated at 2.5 µL extended length
	Racked	FB78010	02-707-440	£18.00	Clear	
	Racked Sterile	FB78012	02-707-441	£20.00	Clear	
0.1 - 20 µL	Bulk	FB78014	02-707-453	£9.00	Clear	 Micropoint tip, Eppendorf style
	Racked	FB78016	02-707-433	£18.00	Clear	
	Racked Sterile	FB78018	02-707-434	£20.00	Clear	
1 - 200 µL	Bulk	FB78020	02-707-450	£9.00	Clear	 Beveled tip
	Racked	FB78022	02-707-418	£18.00	Clear	
	Racked Sterile	FB78024	02-707-419	£20.00	Clear	
	Reload	FB78026	02-707-420	£15.00	Clear	
	Bulk	FB78028	02-707-451	£9.00	Yellow	
	Racked	FB78030	02-707-421	£18.00	Yellow	
	Racked Sterile	FB78032	02-707-422	£20.00	Yellow	
	Reload	FB78034	02-707-423	£15.00	Yellow	
1 - 200 µL	Bulk	FB78036	02-707-446	£10.00	Clear	 Thin Wall Micropoint tip
	Racked	FB78038	02-707-424	£19.00	Clear	
	Racked Sterile	FB78040	02-707-425	£21.00	Clear	
	Reload	FB78042	02-707-426	£15.50	Clear	
	Bulk	FB78044	02-707-452	£10.00	Yellow	
	Racked	FB78046	02-707-427	£19.00	Yellow	
	Racked Sterile	FB78048	02-707-428	£21.00	Yellow	
	Reload	FB78050	02-707-429	£15.50	Yellow	
1 - 200 µL	Bulk	FB78060	02-707-449	£9.00	Clear	 Beveled tip graduated at 10, 50, 100 µL
	Racked	FB78062	02-707-415	£18.00	Clear	
	Racked Sterile	FB78064	02-707-416	£20.00	Clear	
	Reload	FB78066	02-707-417	£15.00	Clear	
	Bulk	FB78052	02-707-448	£9.00	Yellow	
	Racked	FB78054	02-707-412	£18.00	Yellow	
	Racked Sterile	FB78056	02-707-413	£20.00	Yellow	
	Reload	FB78058	02-707-414	£15.00	Yellow	
5 - 300 µL	Bulk	FB78068	02-707-447	£10.50	Clear	 Beveled tip
	Racked	FB78070	02-707-409	£20.50	Clear	
	Racked Sterile	FB78072	02-707-410	£22.50	Clear	
	Reload	FB78074	02-707-411	£17.00	Clear	
100 - 1250 µL	Bulk	FB78076	02-707-400	£12.00	Clear	 Micropoint tip graduated at 100, 200, 500, 1000 µL
	Racked	FB78078	02-707-401	£24.00	Clear	
	Racked Sterile	FB78080	02-707-402	£27.00	Clear	
	Reload	FB78082	02-707-403	£20.00	Clear	
	Bulk	FB78084	02-707-405	£12.00	Blue	
	Racked	FB78086	02-707-406	£24.00	Blue	
	Racked Sterile	FB78088	02-707-407	£27.00	Blue	
	Reload	FB78090	02-707-408	£20.00	Blue	
10 - 20 µL 200 - 300 µL 1250 µL	Empty Rack	FB78092	02-707-445	£12.00		
	Empty Rack	FB78094	02-707-444	£12.00		
	Empty Rack	FB78096	02-707-443	£12.00		

# Features and Benefits


## SureOne Rack Features



- 1 Flush Fit Inserts sit flush with rack, allowing easy multichannel pipetting.
- 2 Recycle Logo molded into all rack components, allows recycling companies to properly recycle the plastic.
- 3 Notched Clasp designed to prop lid open for autoclave venting.
- 4 Sure Stack Lid rim prevents racks from sliding around when stacked and carted around a lab.
- 5 Safe Lock Latch allows for one-handed opening and closing ability and won't pop open when dropped.
- 6 Sure Snap design securely holds insert in place, preventing it from being lifted out while pipetting.
- 7 Flip-Lid can be either a hinged or lift-off lid; hinged format does not prevent use of a multichannel pipet.



## SureOne Insert Features



- 1 Contact Clip prevents unnecessary contact with the tips while clipping in inserts.
- 2 Lot Numbered Inserts for traceability include tip volume and Fisher Cat. No.
- 3 Indented Insert Holes inhibit tip movement, stabilizing them for multichannel work.
- 4 Rigid Inserts are all 96-well multichannel format, offering a robust non-pliable design providing secure seals on every tip.



## SureOne Reload Packaging Features



- 1 Product information displayed on box
- 2 Dust Flap prevents debris from contaminating the product
- 3 Finger Knockouts allow easy access to slide out product
- 4 Canyon Cut Package allows for easy tip insert removal



### Terms and Conditions

Please note that all sales are subject to our standard terms and conditions. Full details can be found on the inside back cover of our catalogues or visit our website [www.fisher.co.uk](http://www.fisher.co.uk)



**Fisher Scientific**

Part of Thermo Fisher Scientific

For Customer Service, call 01509 231166  
 To fax an order, use 01509 231893  
 To order online: [www.fisher.co.uk](http://www.fisher.co.uk)  
 Email: [fsuk.sales@thermofisher.com](mailto:fsuk.sales@thermofisher.com)

Fisher Scientific  
 Bishop Meadow Road, Loughborough  
 Leicestershire, LE11 5RG

©2008 Thermo Fisher Scientific Inc. All rights reserved