

Spoon spatula

The spoon spatula is a double spatula with a 0.5 ml spoon and a 17 mm spoon blade.

The spoon spatula has many different uses in labs – as a sampling spoon, for measuring or transferring small amounts, as a scraper or for mixing and stirring.

Specially designed for disposable use, resistant to chemicals and corrosion-proof.

Available either in a large package with up to 100 pieces or packaged individually and sterilised by gamma rays.

- ▶ PS
- ▶ Production, assembly and packaging according to clean room class 7
- ▶ Complies with EU food requirements and FDA regulations
- ▶ Optionally packaged individually and sterilised by gamma rays

Spoon spatula				
Design	Pieces per pack	Item no.		
Large package	100	5378-0031	N4	
Sterilised, packaged individually	100	5378-1031	N4	



Micro-spatula

Double scoop spatula with rounded ends, content 1.0 and 0.25 ml.

The micro spatula can be used to sample, transfer or process small amounts of chemicals, powders, granulates, pastes, creams or liquids.

The sturdy and convenient centre handle facilitates precise work.

Available either in a large package with up to 100 pieces or packaged individually and sterilised by gamma rays.

- ▶ PS
- ▶ Production, assembly and packaging according to clean room class 7
- ▶ Complies with EU food requirements and FDA regulations
- ▶ Optionally packaged individually and sterilised by gamma rays

Micro-spatula				
Design	Pieces per pack	Item no.		
Large package	100	5378-0032	N4	
Sterilised, packaged individually	100	5378-1032	N4	

