

Bags, polypropylene, biohazard, autoclavable, roll packs

Fisherbrand®



46



Standard temperature bags

- 40µm thickness
- With printed biohazard symbol
- Supplied on convenient rolls of 25 tear off bags per roll
- Suitable for temperatures of up to 121°C

Standard temperature (121°C), HDPE bags. Rolls of 25 bags

Catalogue No	Width, mm	Depth, mm	Capacity, L	No. of rolls
FB55501	250	400	3	40
FB55502	300	500	9	30
FB55503	400	780	21	16
FB55504	500	800	33	20

For autoclaves, please see the Fisher Scientific Laboratory Catalogue or visit our website.



Bags, biohazard

Sterilin



56



Non-sterile sturdy disposal bags with biohazard labelling.

Available in standard temperature, 121°C, or high temperature, 134°C versions. All 'tissue box' style dispensers. Opening along width.

Standard temperature (121°C), polyethylene bags

Catalogue No	Alt. No	Width, mm	Depth, mm	Pack qty
BAJ-440-010Y	509	305	660	200
BAJ-440-020F	509L	405	610	200
BAJ-440-026Q	510L	405	610	500
BAJ-440-040W	511	610	810	200
BAJ-440-051R	510	305	660	500

High temperature (134°C), polypropylene bags

Catalogue No	Alt. No	Width, mm	Depth, mm	Pack qty
BAJ-440-030C	511HT	610	810	200
BAJ-440-046K	510HT	305	660	500
BAJ-440-056H	510LHT	405	610	500
BAJ-440-060Q	509HT	305	660	200
BAJ-440-070N	509LHT	405	610	200

Accessories

Catalogue No	Alt. No	Description	Pack qty
BAJ-440-500X	S23C	Non-sterile metal holder for bag size 305mm x 660mm and 405mm x 610mm	1
BAJ-440-510L	S23B	Non-sterile cardboard holder for bag size 610mm x 810mm	10
BAJ-440-520Y	S23E	Non-sterile metal holder for bag size 610mm x 810mm	1