

Flasks, Erlenmeyer, baffled, Pyrex, glass

PYREX®

NEW

01

For suspension cultures.

- Baffled flask
- Narrow neck
- Borosilicate glass

Catalogue No	Alt. No	Capacity, mL	Pack qty
FHB-276-005H	1134/08	250	5
FHB-276-010X	1134/10	500	5
FHB-276-015E	1134/12	1,000	2
FHB-276-020L	1134/14	2,000	1



Flasks, shaker, borosilicate glass, three baffles

WHEATON
Science International

a

03

Designed for use on either rotary or reciprocating shakers for cell or bacterial cultures.

Special design with three baffles maximises aeration.

Provided with break-resistant lip and marking spot on side of flask. Break-resistant lip accommodates cotton plug or stopper. Wheaton-33® low extractable borosilicate glass (meets ASTM Type I Class A and USP Type I standards).

Catalogue No	Alt. No	Capacity, mL
FHB-450-030T	353255	50
FHB-450-050N	353257	125
FHB-450-070H	353259	250
FHB-450-090B	353262	500
FHB-450-110V	353264	1,000
FHB-450-130P	353266	2,000



For autoclaves, please see the Fisher Scientific Laboratory Catalogue or visit our website.

