

Inoculating loops, sterile



03



Manufactured to the highest quality using cadmium free pigments.


10µL green and 1µL clear loops are manufactured in hard crystal polystyrene for general plating. 10µL white, 5µL red and 1µL yellow are manufactured out of soft HI polystyrene, making them ideal for clinical use as there is no danger of the head breaking off in use.

Catalogue No	Colour	Capacity, µL	Inner pack qty	Pack qty
DIS-105-010M	White	10	20	1,000
DIS-105-020J	Green	10	20	1,000
DIS-105-030G	Red	5	20	1,000
DIS-105-040D	Clear	1	20	1,000
DIS-105-050A	Yellow	1	20	1,000

Individually wrapped

Catalogue No	Colour	Capacity, µL	Pack qty
DIS-105-015C	White	10	1,000
DIS-105-025W	Green	10	1,000
DIS-105-035T	Red	5	1,000
DIS-105-045Q	Clear	1	1,000
DIS-105-055N	Yellow	1	1,000


For burners, please see the Fisher Scientific Laboratory Catalogue, or visit our website.



Inoculating loops, polystyrene, sterile



03



High impact, polystyrene, disposable inoculating loops. 200mm long to reduce the risk of contamination.

Catalogue No	Alt. No	Capacity, µL	Inner pack qty	Pack qty
LPS-131-011B	PL2010/1000	10	20	1,000
LPS-131-021V	PL205/1000	5	20	1,000
LPS-131-031S	PL201/1000	1	20	1,000

Inoculating loops, calibrated



04

Accurately calibrated nichrome wire loops, standard gauge.

Supplied in packs of 25.

Catalogue No	Capacity, µL
LPS-110-010A	10
LPS-110-011V	5
LPS-110-012T	3.3
LPS-110-013R	2.5
LPS-110-014P	2
LPS-110-015N	1

Inoculating loops, heavy gauge



04



Supplied in packs of 25.

Catalogue No	Capacity, µL
LPS-120-010Q	10
LPS-120-011X	5
LPS-120-012M	2.5